cottons

be natural

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Internet: www.cottons.com.au

GB APPLICATOR TAMPONS

1. WASH YOUR HANDS THOROUGHLY BEFORE AND AFTER INSERTING A TAMPON

2. REMOVE THE WRAPPER

Make sure the string is tied securely and hanging outside the applicator. The mechanism of the applicator is similar to that of a syringe. Avoid placing the applicator and tampon on any surfaces after you have unwrapped it.

3. PREPARE TO INSERT

Adopt the position most comfortable for you (eg sit with knees apart, squat or stand with one leg on the toilet, chair or the edge of the bathtub, refer diag. 1) Hold the tampon with the base of the applicator between your thumb and index finger. With the other hand, open the folds of skin at the vaginal opening gently (refer diag. 2) and place the rounded tip of the applicator in the vaginal opening.

4. INSERT THE TAMPON

Breathe out slowly and relax. Then insert the rounded tip of the applicator in, at a slant towards your back, rotating it slightly.

Slide the external insertion tube completely into your vagina, until your fingers touch your body (refer diag. 3a).

Then push the internal insertion tube with your index finger inside until it is inside the first tube. This will position the tampon correctly inside your vagina. Remove both applicator tubes from your body and put them in a wastepaper basket. Make sure the tampon withdrawal string is hanging outside your body. When the tampon is correctly placed you will not feel anything. But if you are unsure, refer to number 5 below.

5. HAVE I INSERTED IT CORRECTLY?

A tampon cannot be felt if it has been inserted correctly, so if you feel any pressure or discomfort the tampon may not have been inserted properly or be far enough inside the vagina. Relax and push it in further or simply wait until the tampon has become saturated so you can easily remove it and start again. Do not insert a tampon if it hurts to do

so and never insert more than one tampon at a time.

6. WHEN TO CHANGE THE TAMPON

By lightly tugging on the withdrawal cord, you will be able to tell if it is time to change – if the tampon does not move, it is not saturated and there is no need to change. Though your needs will depend on your individual flow, you should change your tampon 3-6 times a day (or every 4 to 8 hours) depending on flow.

You certainly shouldn't leave your tampon in for more than 8 hours. We recommend that on each day of your period you use the lowest absorbency tampon needed to suit that day's anticipated menstrual flow.

You should also consider using pads or a sanitary towel overnight. If tampons are worn overnight, insert a fresh tampon before going to sleep and then again first thing in the morning.

7. REMOVING THE TAMPON

Simply pull the string at the same angle at which the tampon was inserted. Check on removal that the tampon is intact. If not, consult your Doctor promptly for advice and assistance. If you ever have trouble finding the string, it can usually be easily reached from a squatting position, using your fingers. Remember to remove the current tampon before inserting another and make sure you remove the last tampon at the end of your period.

8. DISPOSAL OF THE APPLICATOR, TAMPON AND WRAPPER

In the interest of the environment and keeping waterways clean, we recommend you use an alternative means of disposal other than the sewerage system.

9. WHICH SIZE SHOULD I USE?

Some periods are heavy and others light, young girls differ from women after childbirth and a woman's choice of contraception can affect menstruation. So Cottons Applicator Tampons were designed in a range of absorbencies.

COTTONS APPLICATOR TAMPON RANGE AND ABSORBENCIES

REGULAR 

Light to medium flow. 6-9g

SUPER 

Medium to heavy flow. 9-12g

GB TAMPONS

1. WASH YOUR HANDS THOROUGHLY BEFORE AND AFTER INSERTING A TAMPON

2. REMOVE THE WRAPPER

- Remove the coloured tear strip and unwrap the bottom half of the plastic cover. Carefully unwind the withdrawal cord.
- While holding the string firmly, remove the remaining piece of the plastic cover. Try not to handle the main part of the tampon more than is necessary and avoid placing the tampon on any surfaces after you have unwrapped it.

3. PREPARE TO INSERT

Hold the tampon at the cord end. Adopt the position most comfortable for you (eg sit with knees apart, squat, or stand with one leg on the toilet, chair or edge of the bathtub). (Diag 1)

4. INSERT THE TAMPON

(a) Breathe out slowly and relax. With your free hand, spread open the folds of skin at the vaginal opening. (Diag. 2)

(b) Gently push the tampon inside the vaginal opening, the resistance you may feel when you try this is the normal pressure of muscles which will help to hold the tampon in place once you have inserted it. Gently rotating the tampon might help. (Diag. 3b)

(c) Once the tampon is partly inside the vaginal opening, use your middle finger to push the tampon fully inside the vagina. You should aim your finger towards your lower back when guiding the tampon in. Push your finger inside to its full extent, which will place the tampon in the central part of the vagina where there are few nerve endings. Make sure you leave the cord hanging outside your body.

If you are using tampons for the first time, you might need a little practice. It may help to try them for the first time on a heavy day of your period when the vagina will be well lubricated to help insertion. Your finger is just the right length to place a tampon in the correct position inside the vagina. It's just a matter of learning to know your body.

NOTE: Don't worry, a tampon cannot fall out because it is held in place by the muscles at the vaginal opening. It can't get lost either because the opening of the cervix is far too small to allow a tampon to pass through. (Diag. 4)

5. HAVE I INSERTED IT CORRECTLY?

A tampon cannot be felt if it has been inserted correctly, so if you feel any pressure or discomfort the tampon may not have been inserted properly or be far enough inside the vagina. Relax and push it in further or simply wait until the tampon has become saturated so you can easily remove it and start again. Do not insert a tampon if it hurts to do so and never insert more than one tampon at a time.

6. WHEN TO CHANGE THE TAMPON

By lightly tugging on the withdrawal cord, you will be able to tell if it is time to change – if the tampon does not move, it is not saturated and there is no need to change. Though your needs will depend on your individual flow, you should change your tampon 3-6 times a day (or every 4 to 8 hours) depending on flow. You certainly shouldn't leave your tampon in for more than 8 hours. We recommend that on each day of your period you use the lowest absorbency tampon needed to suit that day's anticipated menstrual flow. You should also consider using pads or a sanitary towel overnight. If tampons are worn overnight, insert a fresh tampon before going to sleep and then again first thing in the morning.

7. REMOVING THE TAMPON

Simply pull the string at the same angle at which the tampon was inserted. Check on removal that the tampon is intact. If not, consult your Doctor promptly for advice and assistance. If you ever have trouble finding the string, it can usually be easily reached from a squatting position, using your fingers. Remember to remove the current tampon before inserting another and make sure you remove the last tampon at the end of your period.

8. DISPOSAL OF THE TAMPON AND WRAPPER

In the interest of the environment and keeping waterways clean, we recommend you use an alternative means of disposal other than the sewerage system.

9. WHICH SIZE SHOULD I USE?

Some periods are heavy and others light, young girls differ from women after childbirth and a woman's choice of contraception can affect menstruation. So Cottons Tampons were designed in a range of absorbencies.

COTTONS TAMPON RANGE AND ABSORBENCIES

REGULAR 

Light to medium flow. 6-9g

SUPER 

Medium to heavy flow. 9-12g

GB ABOUT TSS

Cottons 100% natural cotton tampons are safe, hygienic and comfortable to use. They are also hypo-allergenic. You should note though, that no tampons are completely sterile and neither are your hands or vagina, so you should wash your hands thoroughly before and after inserting a tampon.

However tampons, while containing very small amounts of bacteria normally present in the air, have not been shown to carry the bacteria which cause Toxic Shock Syndrome (TSS).

WHAT IS TSS?

Toxic Shock Syndrome (TSS) is a rare but serious illness that may cause death. It is caused by a toxin (a kind of biological poison) which is produced by a type of bacteria (*Staphylococcus aureus*). These bacteria are found in the nose of about one-third of the population. They may also be found on the skin and occasionally in the vagina without causing harm.

WHO IS AT RISK?

TSS can occur in both males and females of any age, but is more common in young women who use tampons during their periods.

HOW DOES IT OCCUR?

If the toxin is produced in the vagina or a wound, and absorbed from there into the bloodstream, a person who is not resistant to the toxin may become ill.

Most people develop resistance to the toxin (that is why the illness is so rare) and in these people there is no harmful effect.

The symptoms of TSS may develop rapidly. Early recognition and treatment of these symptoms can usually prevent serious illness.

WHAT ARE THE SYMPTOMS OF TSS?

The early symptoms may begin suddenly and are similar to the 'flu'. Early recognition of these symptoms is very important.

- Feeling very ill, headache, a sore throat, muscular pains
- High fever and chills, usually 39°C (102°F) or higher
- Vomiting, diarrhoea, or both
- Fainting, dizziness, weakness or confusion
- Sunburn-like rash

WHAT MUST I DO IF I THINK I HAVE TSS?

If, during your period or shortly after, you have any of the above symptoms, REMOVE YOUR TAMPON AND SEE A DOCTOR IMMEDIATELY. Remember it is very important to tell the doctor that you have been using tampons.

If you have ever had TSS, you should not use tampons until you have discussed the matter with a doctor. You may not have developed resistance to the toxin and could get TSS again.

DO TAMPONS CAUSE TSS?

The simple answer is no. Tampons do not carry the bacteria which cause TSS. However, tampon use has been associated with an increased risk of TSS. Although TSS can occur with the use of tampons of any absorbency, the risk increases with the use of tampons of higher absorbency.

WHERE CAN I GET MORE INFORMATION ABOUT TSS?

More information is available, free of charge from the Toxic Shock Syndrome Information Service.

Visit www.toxicshock.com

AS/NZS 2869:1998

GB COTTONS CUSTOMER CARE

In the unlikely event of having to complain about this product, please retain the outer box and send to:

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