# INVENTUM

SINCE 1968

# PRECISION

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# ENGLISH

Legislation requires that all electrical and electronic equipment must be collected for reuse and recycling. Electrical and electronic equipment marked with the symbol indicating separate collection of such equipment must be returned to a municipal waste collection point.


The scale may not be used by persons with medical implants (e.g. pacemakers). Otherwise their function may be impaired.

## DESCRIPTION OF THE APPLIANCE

1. SAFETY INSTRUCTIONS
  2. INTERESTING FACTS
  3. BEFORE FIRST USE
  4. PERSONAL DATA
  5. WEIGHT ONLY
  6. BODY ANALYSIS
  7. EVALUATION OF RESULTS
  8. TIPS
  9. INCORRECT MEASUREMENTS
  10. CLEANING & MAINTENANCE
- GUARANTEE & AFTER SALES SERVICE

PAGE 19  
PAGE 19  
PAGE 20  
PAGE 20  
PAGE 20  
PAGE 20  
PAGE 21  
PAGE 22  
PAGE 22  
PAGE 22  
PAGE 25

1. Glazen weegplateau
2. Meet electrodes
3. Display
4. Insteltoetsen ▲ SET ▼  
Batterijklepje (achterkant)

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
1. Glass surface
2. Measuring electrodes
3. Display
4. Keys for setting ▲ SET ▼  
Battery compartment (bottom)

## 1 | SAFETY INSTRUCTIONS

- **Please read the instructions carefully and completely before operating this appliance and retain for future use.**
- The appliance is intended solely for domestic use within the home.
- Always use batteries which corresponds with the voltage printed on the bottom of the appliance.
- **Attention:** Do not use during pregnancy. Inaccurate measurements may be caused by the amniotic fluid.
- **Attention:** The scale may not be used by persons with medical implants (e.g. pacemakers). Otherwise their function may be impaired.
- **Caution!** Do not stand on the outer edge to one side of the scale otherwise it may tip!
- An electrical appliance is not a toy and should therefore always be placed outside of the reach of children. Children are not always able to understand potential risks. Teach children a responsible handling of electrical appliances.
- Children should be supervised to ensure that they do not play with the appliance.
- This appliance is not intended for use by persons (including children) with reduced physical, sensory or mental capabilities, or lack of experience and knowledge, unless they have been given supervision or instruction concerning use of the appliance by a person responsible for their safety.
- If possible, always weigh yourself at the same time of day (preferably in the morning), after using the toilet, on an empty stomach and unclothed, in order to achieve comparable results.
- Stand on the scale without moving and with equal weight distribution on both feet!
- Minor deviations are normal. When you drink 1 glass of water for example a difference of 300 grams can occur.
- Protect the scale against hard knocks, moisture, dust,

- chemicals, major temperature fluctuations and heat sources which are too close (stove, heating radiators).
- Always ensure that the scale is placed on a flat, stable surface and only use it indoors.
- If you decide never to use the appliance again, we advise you to dispose of it in special containers or at a municipal waste collection point.
- **Never immerse the scale in water or other liquids.**
- The scale can accept a maximum load of 150 kg (330 lb, 24St). The weight measurement is displayed in increments of 100 g (0.2lb, ¼St). The measuring results for the body fat, body water and muscle percentages are displayed in 0.1% increments.
- On delivery the scale is set to the units "cm" and "kg". For changing the measuring unit, please read the information in chapter BEFORE FIRST USE.
- When you want to store the scale for a longer period, we advise to keep it in the original carton.
- Repairs may only be carried out by Inventum customer service or authorized dealers. Before submitting any complaint, first check the batteries and replace them if necessary.
- If a battery has leaked, put on protective gloves and clean the battery compartment with a dry cloth.
- The device is only for private use and is not intended for medical or commercial purposes.
- Please note that technically related measuring tolerances are possible, as the scale is not calibrated for professional medical applications.
- If the appliance is used for other purposes than intended, or is it not handled in accordance with the instruction manual, the full responsibility for any consequences will rest with the user. Any damages to the product or other things are not covered by the warranty.

## 2 | INTERESTING FACTS

### *The measuring principle of the diagnostic scale*

This scale determines the body weight, bodyfat, body water and the percentage of muscle according to the principle of Bioelectric Impedance Analysis (BIA). With this process, an imperceptible, completely safe and harmless current enables the determination of body tissue proportions within seconds. This measurement of electrical resistance (impedance) and the inclusion of constants and/or individual values (age, height, gender) can be used to determine the percentage of fat, water and muscle in the body. Muscle tissue and water have good electrical conductivity, as the fat cells hardly conduct the current due to their very high resistance.

### *Weighing correctly*

If possible, always weighs yourself at the same time of day (ideally in the morning), after going to the toilet, on an empty stomach and without clothes to achieve comparable results. Minor deviations are normal. When you drink 1 glass of water for example a difference of 300 grams can occur. Place the scale on a firm, flat surface; a firm floor covering is a basic requirement for correct weight measurement.

Your scale is an accurate instrument. The human body is the hardest thing to weigh, because it is always in motion. Therefore a 99% accuracy of the body weight is the highest possible. Stand on the scale without moving and with equal weight distribution on both feet! Do not underestimate the weight of clothing and when you have just eaten or drunken anything. Only the long-term trend counts.

## 3 | BEFORE FIRST USE

Before you use the appliance for the first time, please do as follows: carefully unpack the appliance and remove all the packaging material. Keep the material (plastic bags and cardboard) out of reach of children. Check the appliance after unpacking for any damage, possible from transportation. Place the scale on a firm, flat surface. Remove the batteries from the protective packaging and insert with the correct polarity. See the diagram inside the battery compartment.

Press the **UNIT** button on the back of the scale to select the correct measuring unit (kg/st/lb).

kg = kilograms

lb = pounds

st = stone

Place the scale on a flat and stable surface and wait for the scale to turn itself off.

## 4 | PERSONAL DATA

You must enter personal user parameters in order to measure your body fat percentage and other physical data. The scale has 8 memory positions in which you and other members of your family can save and retrieve personal settings.

Now you can enter the following setting:

Memory position (P)	1 to 8
Body size	80 to 220 cm
Age	6 to 100 years
Gender	male (♂), female (♀)
Sportsman or -woman	♣

- Press the **SET** key when the display is switched off to change the setting. The first memory position **P1** flashes on the display. To select the correct user, press the **▲** or **▼** keys and press the **SET** key again to confirm.
- After you have entered the correct user, the symbol for gender will blink. To select the gender (male, female, sportsman, sportswoman), press the **▲** or **▼** keys and press the **SET** key again to confirm.
- After you have entered the gender, the indication for body size (height) will blink (preset 165 cm). To set your own height press the **▲** or **▼** keys and press the **SET** key again to confirm.
- After you have entered your body size, the indication for age will blink (preset 30). To set the correct age (male, female, sportsman, sportswoman), press the **▲** or **▼** keys and press the **SET** key again to confirm.

After all parameters have been entered, you can now measure your weight, body fat and the other data. The display will show **0,0 kg**.

Stand on the scale barefoot and let the scale perform a measuring. When you step on the scale the next time the scale will recognize you. When you step off the scale, the scale will switch off after a few seconds.

*TIP: To reset a memory position for another user perform the same steps as described above. All positions are preset.*

## 5 | WEIGHT ONLY

The scale has a step-on function. Step barefoot onto the scale. Stand still with your weight distributed evenly on both legs. The scale will start measuring your weight immediately and a short time later will display your weight. When you step off the scale it will turn off after a few seconds.

## 6 | BODY ANALYSIS

When you stay on the scale after it displayed your weight, the scale will automatically perform a body analysis. The scale can identify users according to the weight.

*Important: When using the analysis function, make sure there is no contact between your feet, legs, calves or thighs. Otherwise the measurement will not be performed correctly.*

After the weight is stable, the "0" signal moves from left to right for two times, and then shows the data.

The scale is measuring the bodyfat, -water, etc.

The scale can identify a person by their weight. If it detects more users with a similar weight, it will display the relative user positions. Use the ▲ or ▼ keys to choose the right user. Stay on the scale to complete the measurement.

The scale will now display the results:

bodyfat

bodywater

muscle

and your bone mass will be shown in kilograms.

To give you the opportunity to see all the results, the scale will show them three times before it turns off:

1. weight
2. bodyfat
3. bodywater
4. bodymuscle
5. bone mass

If you don't choose a user or when the scale can't identify you, it will only show your weight and then turn off.

## 7 | EVALUATION OF RESULTS

The percentage of bodyfat represents the in terms of percentage part of the complete body weight. There is no real advise for this value, because it depends on gender and age.

Bodyfat is very important for your body, to ensure good functioning. It protects organs, hinges, regulates the body temperature, stores vitamins and acts as a reserve for the body. Females need a bodyfat percentage between 10% and 13%. Males need a bodyfat percentage between 2% and 5%, with exception of sports people.

The following **body fat** percentages are guide values (contact your physician for further information).

Man

Age	-    ○    +    ++			
	low	normal	high	very high
10-16	< 10	10 - 18	19 - 23	> 23
17-39	< 12	12 - 20	21 - 25	> 25
40-55	< 13	13 - 21	22 - 26	> 26
56-100	< 14	14 - 22	23 - 27	> 27

Woman

Age	-    ○    +    ++			
	low	normal	high	very high
10-16	< 18	18 - 28	29 - 35	> 35
17-39	< 20	20 - 32	33 - 38	> 38
40-55	< 23	23 - 35	36 - 41	> 41
56-100	< 24	24 - 36	37 - 42	> 42

A lower value is often found in athletes. Depending on the type of sports, training intensity and physical constitution, values may result which are below the recommended values states. It should, however, be noted that there could be a danger to health in the case of extremely low values.

A large part of the human body contains water. This **bodywater** can be found through the entire body. About 60% of the bodymass consists of water with men. With women this is about 55% (caused by a higher bodyfat percentage).

The **body water** percentage is normally within the following ranges:

Man

Age	poor	good	very good
10 - 15	< 58	58.1 - 72	> 72
16 - 30	< 53	53.1 - 67	> 67
31 - 60	< 47	47.1 - 61	> 61
61 - 100	< 42	42.1 - 56	> 56

Woman

Age	poor	good	very good
10 - 15	< 57	57.1 - 67	> 67
16 - 30	< 47	47.1 - 57	> 57
31 - 60	< 42	42.1 - 52	> 52
61 - 100	< 37	37.1 - 47	> 47

Body fat contains relatively little water. Therefore persons with a high body fat percentage have body water percentages below the recommended values. With endurance athletes, however, the recommended values could be exceeded due to low fat percentages and high muscle percentage.

This scale is unsuitable for measuring body water in order to draw medical conclusions concerning age-related water retention, for example. If necessary ask your physician.

The **muscle percentage** is normally within the following ranges:

Age	MAN healthy	WOMAN healthy
10-99	> 40	> 34

Like the rest of our body, our bones are subject to the natural development, degeneration and ageing processes. Bone mass increases rapidly in childhood and reaches its maximum between 30 and 40 years of age. Bone mass reduces slightly with increasing age. You can reduce this degeneration somewhat with healthy nutrition (particularly calcium and vitamin D) and regular exercise. With appropriate muscle building, you can also strengthen your bone structure. Note that this scale will not show you the calcium content of your bones, but will measure the weight of all bone constituents (organic substances, inorganic substances and water).

The bone mass is normally within the following ranges:

	MAN			WOMAN		
WEIGHT	< 60 kg	60-75 kg	> 75 kg	< 45 kg	45-60 kg	> 60 kg
BONE MASS	2,5 kg	2,9 kg	3,2 kg	1,8 kg	2,2 kg	2,5 kg

**Attention:** Please do not confuse bone mass with bone density. Bone density can be determined only by means of a medical examination (e.g. computer tomography, ultrasound). It is therefore not possible to draw conclusions concerning changes to the bones and bone hardness (e.g. osteoporosis) using this scale.

Little influence can be exerted on bone mass, but it will vary slightly within the influencing factors (weight, height, age, gender).

## 8 | TIPS

Remember that only long-term trends are important. Short-term fluctuations in weight over a few days are usually the result of a loss of fluid.

Important when measuring body fat/body water/muscle percentage:

- If possible, always weighs yourself at the same time of day (ideally in the morning), after going to the toilet, on an empty stomach and without clothes to achieve comparable results.
- The measurement may only be carried out while barefoot and it is helpful if the soles of your feet are slightly damp. Completely dry soles can result in unsatisfactory results, as they have insufficient conductivity.
- Stand still during the measurement.
- Wait several hours after unusually strenuous activity.
- Wait approx. 15 minutes after getting out of bed so that the water in your body can be distributed.
- Remember that only long-term trends are important. Short-term changes in weight within a few days are usually caused by a loss of fluid. Body water plays an important role towards our general wellbeing.

The measurement is not reliable for:

- Children under approx. 10 years of age;
- Professional sporters and bodybuilders;
- Pregnant women;
- Persons with fever, undergoing dialysis, with symptoms of edema or osteoporosis;
- Persons taking cardiovascular medication.
- Persons taking vasodilating or vasoconstricting medications;
- Persons with substantial anatomical deviations in the legs relative to their total height (leg length considerably shorter or longer than usual).

## 9 | INCORRECT MEASUREMENTS

If the scale detects an error during measurement, this is indicated with ----.

Possible causes of errors	Causes of error
The electrical resistance between the electrodes and the soles of your feet is too high (e.g. with heavily callused skin).	Repeat weighing barefoot. Slightly moisten the soles of your feet if necessary. Remove the calluses on the soles of your feet if necessary.
The fat percentage is outside the measurable range (less than 3% or greater than 50%).	Repeat weighing barefoot. Slightly moisten the soles of your feet if necessary.
The water percentage is outside the measurable range (less than 25% or greater than 75%).	Repeat weighing barefoot. Slightly moisten the soles of your feet if necessary.
The muscle and bone percentages are outside the measurable range - age and gender related (less than 25% or greater than 75%).	Repeat weighing barefoot. Slightly moisten the soles of your feet if necessary.
The maximum load-bearing capacity of 150 kg (330 lbs) was exceeded ----.	Only weigh the maximum permissible weight.

## 10 | CLEANING & MAINTENANCE

The scale should be cleaned occasionally.

Dirt and dust has a tendency to collect on the rubber electrodes. Always stand on the scale barefoot or wearing socks as the rubber surfaces will quickly become dirty if you wear shoes. Use a damp cloth to clean the scale. Do not use abrasive detergents and never immerse the scale in water.

- **Do not use aggressive solvents or cleaning agents!**
- **Never immerse the scale in water or any other liquids.**
- **Handle the scale with care, protect the scale from impact with hard object and do not drop it.**

Your scale is fitted with a change battery indicator. If the scale is operating with batteries which are too weak, **Lo** appears on the display and the scale will switch off automatically. In this case, the batteries must be exchanged (2 x CR2032, 3V). If necessary please use a small screwdriver to remove the batteries.

Batteries can contain toxins that are harmful to health and the environment. The used, completely drained standard and rechargeable batteries must be disposed of in specially marked collection containers, at toxic waste collection points or electrical product retailers. You are legally obligated to dispose of the batteries.

**Note:** Batteries containing pollutants are marked with the following symbols: Pb = Battery contains lead, Cd = Battery contains cadmium, Hg = Battery contains mercury.

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## GUARANTEE & AFTER SALES SERVICE

- This product is guaranteed for a period of 2 years from the date of purchase against mechanical and electrical defect.
- Guarantee can only be given if proof of purchase (e.g. till receipt) can be shown. Therefore keep it in a safe place. The guarantee period starts at date of purchase.
- Guarantee is not valid in case of:
  - normal wear;
  - damage through misuse or when the instructions were not followed correctly;
  - lack of maintenance;
  - not following the operation- and safety instructions;
  - incompetent montage or repairs made by yourself or others;
  - using accessories and spare parts which are not suitable for the appliance; for assessment by our technical department.
- The guarantee is not valid for:
  - normal wear;
  - batteries;
  - (interchangeable) signal- and controllamps for assessment by our technical department.

Return to your supplier or send the appliance, including a copy proof of purchase and description of the complaint, stamped to:

Inventum Group BV  
Generatorstraat 17  
3903 LH VEENENDAAL  
The Netherlands  
Website: [www.inventum.eu](http://www.inventum.eu)

IES Import Electronic Service GmbH  
Lembekstrasse 36  
22529 Hamburg  
GERMANY  
Website: [www.ihrservice.de](http://www.ihrservice.de)

- Pack it carefully (preferably in the original carton). Ensure unit is clean.
- Enclose your name and address and quote model number on all correspondence.

When the guarantee is no longer valid or when the period of guarantee has passed, we will charge any costs that flow out of an assessment and/or repairing.

On request we will supply a statement before repairing. Without requesting a statement, we assume that you, by sending the appliance, automatically order to repair the appliance.

A guarantee period of two months is valid for repairs done by our technical department and the used spare parts, counting from the repair date. When a malfunction occurs within this period which directly flows out of the repairs and the used spare parts, the consumer has the right of repairs without any costs.

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GL-International GmbH