SANITAS

SBG 43

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ENGLISH

Contents

| | |
|--|----|
| 1. Getting to know your instrument | 15 |
| 2. Signs and symbols | 15 |
| 3. Safety notes | 15 |
| 4. Unit description..... | 17 |
| 5. Initial use | 17 |
| 6. Setting up the control panel | 17 |
| 7. Information..... | 18 |
| 8. Settings..... | 18 |
| 9. Operation | 19 |
| 10. Evaluation of results..... | 21 |
| 11. Replacing batteries | 22 |
| 12. Cleaning and care the unit..... | 23 |
| 13. Disposal | 23 |
| 14. What if there are problems? | 24 |
| 15. Technical specifications | 25 |

Included in delivery

- Diagnostic scale
- Removable control unit with display
- Wall mount for attachment
- 4 x 1.5 V AAA batteries (scale)
- 2 x CR2032 lithium batteries (control unit)
- These instructions for use

Please read these instructions for use carefully and keep them for later use, be sure to make them accessible to other users and observe the information they contain. This diagnostic scale is designed to make a significant contribution to your health.

1. Getting to know your instrument

Function of the unit

This digital diagnostic scale is used to weigh you and to provide diagnostic data about your personal fitness. The scale has a display that receives the measured values by infrared transmission. Weight is displayed in 100 g graduations. The scale is designed for private use.

The scale has the following diagnostic functions that can be used by up to 10 persons:

- Body weight measurement
- Measurement of body fat percentage
- Body water percentage
- Muscle percentage
- Bone mass
- Basic and active metabolic rate

The memory function allows you to access the two most recent measurements.

This scale also has the following functions:

- Interpretation of body fat values
- Switching between kilograms "kg" and pounds "lb"
- Automatic shutoff function
- Battery change indicator for weak batteries

2. Signs and symbols

The following symbols appear in these instructions.

Warning

Warning instruction indicating a risk of injury or damage to health.

Important

Safety note indicating possible damage to the unit/accessory.

Note

Note on important information.

3. Safety notes

Please read these instructions for use carefully and keep them for later use, be sure to make them accessible to other users and observe the information they contain.

WARNING

• **The scales must not be used by persons with medical implants (e.g. heart pacemakers) Otherwise their function could be impaired.**

- Do not use during pregnancy.
- Do not stand on the outer edge to one side of the scale otherwise it may tip!
- Keep children away from packaging materials (risk of suffocation).
- Caution! Do not step onto the scale with wet feet or when the scale's surface is damp - danger of slipping!

Notes on handling batteries

- Swallowing batteries can be extremely dangerous. Keep the batteries and scale out of the reach of small children. Should a child swallow a battery, seek medical assistance immediately.
- Replace weak batteries before they discharge completely.
- Always replace all batteries at the same time and use batteries of the same type.
- Batteries must not be recharged, taken apart, thrown into an open fire or short circuited.
- Leaking batteries may damage the device. If you do not intend to use the device for longer periods, remove the batteries from the battery compartment.
- If a battery has leaked, put on protective gloves and clean the battery compartment with a dry cloth.

- Batteries can contain toxins that are harmful to health and the environment. Always dispose of batteries in accordance with applicable legal regulations. Do not dispose of batteries with the normal household waste.

General notes

- The unit is for personal use only and is not intended for medical or commercial applications.
- Caution! Do not step onto the scale with wet feet or when the scale's surface is damp - danger of slipping!
- Note that technical tolerances mean that results may vary, because the scale is not calibrated for professional medical use.
- The scale can accept a maximum load of 150 kg (330 lb). The weight measurement is displayed in increments of 100 g (0.2 lb).
- The measuring results of the body fat, body water and muscle percentage are shown in 0.1% increments.
- On delivery, the scale is set to „cm“ and „kg“. At the back of the scales is a switch to change between „pounds“ and „stones“.
- Place the scale on a firm level floor; a firm floor covering is required for correct measurement.
- Protect the unit against hard knocks, moisture, dust, chemicals, major temperature fluctuations and heat sources which are too close (stoves, heating radiators).
- Repairs may only be carried out by customer service or authorised dealers. Before submitting any complaint, first check the batteries and replace them if necessary.
- Should you have any questions concerning the use of our devices, please contact your dealer or customer service.
- All scales comply with EC Directive 2004/108/EC. Should you have any questions concerning the use of our devices, please contact your dealer or customer service.

Storage and maintenance

The precision of the measured values and service life of the unit depend on careful use:

IMPORTANT

- The unit should be cleaned occasionally. Do not use abrasive detergents and never immerse the unit in water.
- Ensure that no liquid gets on the scale. Never dip the scale into water. Never wash it under running water.
- Do not place any objects on the scale when it is not being used.
- Protect the scales and the control panel from impact, moisture, dust, chemicals, major temperature fluctuations and make sure it is not too close to heat sources (ovens, radiators). Do not press the buttons using excessive force or with pointed objects.
- Do not expose the scale and the control panel to high temperatures or strong electromagnetic fields (e.g. mobile phones)

Batteries

WARNING

- Swallowing batteries can be extremely dangerous. Keep the batteries out of the reach of small children. If a battery is swallowed, get immediate medical help.
- Do not throw batteries into fire. Explosion hazard!

IMPORTANT

- Batteries must not be recharged or reactivated by other means, taken apart or short-circuited.
- Remove batteries if the unit will not be used for an extended period of time. This prevents damage that may occur from leaks.

Repairs


IMPORTANT

- Never open or attempt to repair the unit yourself, as otherwise proper function is no longer guaranteed. Failure to observe this regulation shall void the warranty.
- If you need to have the unit repaired, contact customer service or an authorised dealer.

4. Unit description

Overview

1. Control panel with display
2. SET/ON button
3. Down key
4. Up key
5. Receiver
6. Scale
7. Transmitter
8. LED
9. Wall mount


5. Initial use

Inserting batteries

If present, pull the battery insulating strip off the battery compartment cover or remove the battery's protective film and insert the battery in the scale and the control panel with the correct polarity. Note the graphic in the battery compartment.

Control unit

2 x CR2032 3 V lithium batteries


Scale

4 x 1.5 V AAA batteries

If the scale fails to operate, remove the batteries completely and reinsert them.

For instructions on changing the batteries, see p. 23.

Changing the weight unit

On delivery the scale is set to the unit „kg“. A switch located on the back of the control panel enables you to change to „pounds“ (lb).

Setting up the scale

Set up the scale on a level and solid surface. A solid surface is a prerequisite for a correct measurement.

6. Setting up the control panel


Positioning the control panel

The control panel can be held in the hand, attached to the wall mount or be left on the scale.

Wall mounting

The wall mount can be attached to the wall with the double-sided adhesive strip. The surface of the wall must be clean, dry, plane and grease-free.

- Remove the foil from the adhesive strip on the back of the wall mount
- Press the wall mount firmly against the wall.
- Make sure that the wall mount is straight.
- Note that the control panel and the scale must also be aligned to each other.

7. Information

The measuring principle

This scale operates according to the B.I.A. principle (bioelectric impedance analysis). This enables the measurement of physical relationships within seconds by means of an undetectable, completely harmless electric current. The body fat percentage and other physical relationships in the body can be determined by measuring the electrical resistance (impedance) and calculating constants and individual parameters such as age, height, gender and degree of activity.

Muscle tissue and water have good electrical conductivity, and therefore lower resistance. Bones and fatty tissue, on the other hand, have low conductivity as fat cells and bones conduct hardly any current as a result of their very high resistance.

Please be aware that values obtained from the diagnostic scale represent only an approximation of actual analytical medical data. Only a specialist physician can accurately determine body fat, body water, muscle percentage and bone structure using medical procedures (e.g. computed tomography).

General tips

- If possible, always weigh yourself at the same time of day (preferably in the morning), after using the toilet, on an empty stomach and unclothed, in order to achieve comparable results.
- Important for the measurement: only measure body fat when barefoot and with the soles of the feet slightly moist.
Completely dry soles can result in unsatisfactory measurements due to inadequate conductivity.
- Stand still and upright during the measurement.
- Wait a few hours after unaccustomed physical activity.
- Wait approx. 15 minutes after getting up in the morning to allow the water in the body to distribute.
- Remember that only long-term trends are important. Short-term changes in weight within a few days are usually caused by a loss of fluid. Body water plays an important role towards our general wellbeing.

Limitations

When measuring body fat and other values, deviating and implausible results may occur in:

- children under approx. 10 years,
- competitive athletes and body builders,
- pregnant women,
- persons with fever, undergoing dialysis treatment or with symptoms of edema or osteoporosis,
- persons taking cardiovascular medication (affecting the heart and vascular system),
- persons taking vascular dilating or vascular constricting medication,
- persons with considerable anatomic deviations of the legs in relation to total body size (length of the legs considerably shortened or lengthened).

8. Settings

Enter your personal information before using the scale.

On delivery, the scale is set to „cm“ and „kg“. At the back of the scales is a switch to change to „pounds“.

Setting user data

You must enter personal user parameters in order to measure your body fat percentage and other physical data.

The scale has 10 user memory locations which allows, for example, you and your family members to save personal settings and retrieve them later.

- Switch the control panel on by pressing the „[SET/ON]“ key.
- Then wait until the display “ - - - - “ appears.
- Start setting with the „[SET/ON]“ key. The first memory space flashes in the display.
- Select the required memory space with the ▼/▲ key and confirm with the Set key
- Now you can adjust the following settings:

| User data | Parameters |
|--------------------|----------------------------------|
| Memory position | 1 to 10 |
| Body height | 100 to 220 cm (3' 03" to 7' 03") |
| Age | 10 to 100 Jahre |
| Sex | man (♂), woman (♀) |
| Degree of activity | 1 to 5 |

- Modifying values: Press ▲ or ▼ or press and hold for more rapid progress.
- Confirming entries: Press the „[SET/ON]“ button.
- Once the values have been displayed, “ - - - - “ appears in the display.
- Audibly, briefly and firmly tap the scale with your foot to activate the vibration sensor. The LED on the scale first shows red and then switches to green after about 5 seconds. At the same time, the display switches to “0.0”.
- Now the scale is ready to measure your weight. Now step onto the scale. Stand still on the scale with your weight distributed evenly between both legs. The scale immediately begins to measure your weight.

If no further action is performed, the scale and the control panel switch off.

Degrees of activity

Selection of the degree of activity must refer to the medium and long term.




| Degree of activity | Physical activity |
|--------------------|---|
| 1 | None. |
| 2 | Low: A small amount of light physical effort (e.g. short walks, light garden work, gymnastic exercises). |
| 3 | Medium: Physical effort for 30 minutes at least 2 to 4 times a week. |
| 4 | High: Physical effort for 30 minutes at least 4 to 6 times a week. |
| 5 | Very high: Intensive physical effort, intensive training or hard physical work for at least one hour daily. |

All parameters must be entered before the body fat and other values can be entered.

9. Operation

Place the scale on a firm level floor (no carpet); a firm floor covering is required for correct measurement.

Make sure that the control panel and the scale are in line of sight to each other. Do not cover the receiver and transmitter. The maximum distance between control panel and scale is 2 m.

Only weight measurement

Press the „[SET/ON]“ button on the display module. The entire display appears as a self-test until “- - - -” is displayed (Fig. 1).

Audibly, briefly and firmly tap the scale with your foot to activate the vibration sensor. The LED first shows red and then switches to green after about 5 seconds.

If the control panel is receiving the signal from the scale it shows „0.0“ (Fig. 2). Now the scale is ready to measure your weight. Now step onto the scale. Stand still on the scale with your weight distributed evenly between both legs. The scale immediately begins to measure your weight.

Fig. 1

Fig. 2

Weighing, carrying out diagnostics

- Press the „[SET/ON]“ button on the control panel. The complete LCD display appears, then “- - - -” is displayed.
- Now press the button ▲.
- Select the memory space where your personal data are stored by repeatedly pressing the ▲ or ▼ keys and confirm with „[SET/ON]“. These data are then displayed consecutively until “- - - -” appears (Fig. 1).
- Audibly, briefly and firmly tap the scale with your foot to activate the vibration sensor. The LED first shows red and then switches to green after about 5 seconds. If the control panel is receiving the signal from the scale it shows „0.0“ (Fig. 2).
- Step onto the scale barefoot and make sure that you are standing still with your weight evenly distributed between both legs on the stainless steel electrodes.

Fig. 3

Fig. 4

(i) Note: Your feet, legs, calves and thighs must not touch each other. Otherwise the measurement cannot be performed correctly. The scale immediately begins to measure your weight.

The scales begin weighing immediately. First, the weight is displayed. Whilst the further parameters are being measured an „arrow on the left edge of the display“ appears (fig. 3). Your measured weight is displayed shortly after this (Fig 4).

Measure weight, call up stored values, perform diagnosis

- Press the „[SET/ON]“ button on the control panel. The complete LCD display appears, then “- - - -” is displayed.
- Now press the button. ▼.
- Select the memory space where your personal data are stored by repeatedly pressing the ▼ or ▲ keys and confirm with „[SET/ON]“.
- Firstly your stored data are displayed consecutively, followed by the last measured values.

Firstly the following display appears.

1. Weight in kg
2. Body fat in % with interpretation of BF
3. Body water in %
4. Muscle percentage in %
5. Bone mass in kg
6. Active metabolic rate in kcal

Switching off the scale

The scale and the control panel switch off automatically after 10 minutes. While you adjust the settings it takes 90 seconds before the control panel automatically switches off.

Hold the „[SET/ON]“ button depressed for 3 seconds to switch off the display module manually.

Make sure that you always switch on the scale first and wait for the display to show „0.0“ (Fig. 2) before stepping on it.

10. Evaluation of results

Body fat percentage

Body fat content is evaluated with the arrow on the left edge of the display. The following body fat percentages are guide values (contact your physician for further information).

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A lower value is often found in athletes. Depending on the type of sports, training intensity and physical constitution, values may result which are below the recommended values stated. It should, however, be noted that there could be a danger to health in the case of extremely low values.

Body water percentage

The body water percentage is normally within the following ranges:

Man

| Age | poor | good | very good |
|--------|------|---------|-----------|
| 10-100 | <50 | 50-65 % | >65 |

Woman

| Age | poor | good | very good |
|--------|------|---------|-----------|
| 10-100 | <45 | 45-60 % | >60 |

Body fat contains relatively little water. Therefore persons with a high body fat percentage have body water percentages below the recommended values. With endurance athletes, however, the recommended values could be exceeded due to low fat percentages and high muscle percentage.

This scale is unsuitable for measuring body water in order to draw medical conclusions concerning age-related water retention, for example. If necessary ask your physician. Basically, a high body water percentage should be the aim.

Muscle percentage

The muscle percentage is normally within the following ranges:

Man

| Age | low | normal | high |
|--------|-----|---------|-------|
| 10-14 | <44 | 44-57 % | >57 % |
| 15-19 | <43 | 43-56 % | >56 % |
| 20-29 | <42 | 42-54 % | >54 % |
| 30-39 | <41 | 41-52 % | >52 % |
| 40-49 | <40 | 40-50 % | >50 % |
| 50-59 | <39 | 39-48 % | >48 % |
| 60-69 | <38 | 38-47 % | >47 % |
| 70-100 | <37 | 37-46 % | >46 % |

Woman

| Age | low | normal | high |
|--------|-----|---------|-------|
| 10-14 | <36 | 36-43 % | >43 % |
| 15-19 | <35 | 35-41 % | >41 % |
| 20-29 | <34 | 34-39 % | >39 % |
| 30-39 | <33 | 33-38 % | >38 % |
| 40-49 | <31 | 31-36 % | >36 % |
| 50-59 | <29 | 29-34 % | >34 % |
| 60-69 | <28 | 28-33 % | >33 % |
| 70-100 | <27 | 27-32 % | >32 % |

Bone mass

Like the rest of our body, our bones are subject to the natural development, degeneration and ageing processes. Bone mass increases rapidly in childhood and reaches its maximum between 30 and 40 years of age. Bone mass reduces slightly with increasing age. You can reduce this degeneration somewhat with healthy nutrition (particularly calcium and vitamin D) and regular exercise.

With appropriate muscle building, you can also strengthen your bone structure. Note that this scale will not show you the calcium content of your bones, but will measure the weight of all bone constituents (organic substances, inorganic substances and water). Little influence can be exerted on bone mass, but it will vary slightly within the influencing factors (weight, height, age, gender). No recognised guidelines or recommendations relating to bone mass measurement.

ATTENTION:

Please do not confuse bone mass with bone density.

Bone density can be determined only by means of a medical examination (e.g. computer tomography, ultrasound). It is therefore not possible to draw conclusions concerning changes to the bones and bone hardness (e.g. osteoporosis) using this scale.

AMR

The active metabolic rate (AMR) is the amount of energy required daily by the body in its active state. The energy consumption of a human being rises with increasing physical activity and is measured on the diagnostic scale in relation to the degree of activity entered (1-5).

To maintain your existing weight, the amount of energy used must be reintroduced into the body in the form of food and drink. If less energy is introduced than is used over a longer period of time, your body will obtain the difference largely from the amount of fat stored and your weight will decrease. If, on the other hand, over a longer period of time more energy is introduced than the total active metabolic rate (AMR) calculated, your body will be unable to burn off the excess energy, and the excess will be stored in the body as fat and your weight will increase.

Results in relation to time

Remember that only long-term trends are important. Short-term fluctuations in weight over a few days are usually the result of a loss of fluid.

The interpretation of the results will depend on changes in your overall weight and body fat, body water and muscle percentages, as well as on the period during which these changes take place. Sudden changes within days must be distinguished from medium term changes (over weeks) and long term changes (months).

A basic rule is that short term changes in weight almost exclusively represent changes in water content, whereas medium and long term changes may also involve the fat and muscle percentages.

- If your weight reduces over the short term, but your body fat percentage increases or remains the same, you have merely lost water – e.g. after a training session, sauna session or a diet restricted only to rapid weight loss.
- If your weight increases over the medium term and the body fat percentage falls or stays the same, then you could have built up valuable muscle mass.
- If your weight and body fat percentage fall simultaneously then your diet is working – you are losing fat mass.
- Ideally you should support your diet with physical activity, fitness or power training. This way you can increase your muscle percentage over the medium term.
- Body fat, body water or muscle percentages should not be added (certain elements of muscle tissue also contain body water).

11. Replacing batteries

The control panel is equipped with a „replace battery“ indicator.

The display shows „L0-1“ and „L0-2“ when the batteries are beginning to run down and the scale switches off automatically. The batteries must be replaced as soon as possible. If the display shows „L0-1“, the bat-

teries in the scale must be replaced (4 x 1.5 V AAA). If the display shows „L3-2“, the batteries in the control panel must be replaced (2 x CR2032 lithium batteries).

The batteries in the scale must also be replaced if the LED on the scale does not light after activating the vibration sensor (4 x 1.5 V AAA batteries).

i NOTE:

- When changing batteries, always use batteries of the same type, brand and capacity.
- Do not use rechargeable batteries.
- Use batteries that are free of heavy metals.

12. Cleaning and care of the unit

The unit should be cleaned occasionally.

Clean using a damp cloth, to which you can apply a little detergent if necessary.

IMPORTANT

- Do not use aggressive solvents or cleaning agents!
- Never immerse the unit in water!
- Do not wash the unit in a dishwasher!

13. Disposal

Standard and rechargeable batteries should not be disposed of separately from the household waste. As a consumer, you are legally obliged to return used batteries for proper disposal. You can hand in your used batteries at public collection points in your district or sales outlets where batteries of this type are sold.

i NOTE:


The codes below are printed on batteries containing harmful substances:

- Pb = Battery contains lead,
Cd = Battery contains cadmium,
Hg = Battery contains mercury.

For environmental reasons, do not dispose of the scale with the batteries in the household waste at the end of its useful life. Dispose of the unit at a suitable local collection or recycling point.

Observe the local regulations for material disposal.

Dispose of the unit in accordance with EC Directive 2002/96/EC – WEEE (Waste Electrical and Electronic Equipment).

If you have any questions, please contact the local authorities responsible for waste disposal.

14. What if there are problems?

If the scale detects an error when measuring, the following is displayed.

| Display | Cause | Remedy |
|------------|---|--|
| 0-Ld | The maximum load-bearing capacity of 150 kg was exceeded. | The load must not exceed 150 kg. |
| Lo-1 | The scale batteries are almost flat. | Replace the batteries of scale. |
| Lo-2 | The control panel batteries are almost flat. | Replace the batteries of control panel. |
| Err | Faulty or no transmission between scale platform and control panel. | Place the scale on a flat surface and activate it correctly (wait for the display to show „0.0“) and repeat the measurement (see page 19). |
| No display | The batteries in the control panel are entirely flat. | Replace the batteries (see page 22). |
| | The batteries in the control panel are not correctly installed. | Check that the polarity is correct (Note the graphic in the battery compartment.). |

No measurement possible

| Possible errors | Remedy |
|--|--|
| The scale was not activated before being stepped on. If you step onto the scale before „0.0“ appears in the display, it will not function correctly. | Activate the scale correctly (wait for „0.0“) and repeat the measurement. Switch on the control panel and the scale before the measurement. If necessary, wait 30 seconds and then repeat the measurement. |
| Interference between scale and control panel. | Remove and replace the batteries in the scale and control panel or make sure that the transmitter and receiver are not covered. |
| The control panel is too far from the scale. | Reduce the distance between the scale and control panel. |
| The scale batteries are flat. | Replace the batteries (see page 22). |
| The connection between control panel and scale is interrupted. | Do not cover the receiver and transmitter. |

Incorrect weight measurement

| Possible errors | Remedy |
|----------------------------------|--|
| Scale is on carpet. | Place scale on a flat, solid surface. |
| Scale zero setting is incorrect. | Switch on scale and wait until it switches off automatically (LED off). Then repeat measurement. |

15. Technical specifications

| | |
|------------------------|------------------------------|
| Scale batteries | 4 x 1.5 V AAA |
| Control unit batteries | 2 x CR2032 lithium batteries |
| Measuring range | 1 kg to 150 kg |
| Weight | 100 g |
| Transmission | Infrared transmission |

Subject to technical modifications.

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