DK24636


Our organic cotton is certified according to the Global Organic Textile Standard (GOTS)

Tampon ingredients:
100% certified organic cotton

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
Your vagina is
one of the most
absorbent parts
of your body.
Let's be extra
careful there.

Join us at
www.yoni.care

If you have any questions or feedback
please contact us at: info@yoni.care

V8.H3818

How to use Yoni tampons?

Step 1 GETTING READY

- ◆ Wash your hands.
- ◆ Take a look at your Yoni tampon. Is it totally sealed? Yes, good. (If not, don't use it and take another one.)
- ◆ Find the colored tear strip and pull. Once you've done that, removing both ends of the transparent wrapper is easy.
- ◆ Find the string at the end of the tampon and pull on it until the end slightly flares out. This widens the base of the tampon a little, creating a nice and protected pocket for your finger. This is handy since your finger will be inserting the tampon.

Step 2 INSERTING THE TAMPON

- ◆ Place your index finger in the little pocket. Now with your thumb and middle finger get a firm hold on the tampon. Next stand with one leg up (rest your foot on e.g. the toilet) or sit with your legs open.
- ◆ Take a deep breath, exhale and relax. With your free hand, gently hold open the skin around the vaginal opening. Now take the Yoni tampon and with the full length of your index finger, ease it into your vagina, aiming toward your lower back.
- ◆ When the tampon is in the right place, you won't feel it. If you do, push it a little further in. Leave the string hanging outside of your body.
- ◆ Wash your hands and you're done!

Step 3 REMOVING THE TAMPON

- ◆ You should regularly change your tampon. This can vary between every 4 to max. 8 hours. If you start to feel your tampon it is a sign that it's time to change it.
- ◆ First wash your hands.
- ◆ Take a deep breath, exhale and relax, gently but firmly tug at the string and it's out! This should be easy to do, otherwise you may be changing your tampon too often or using a too high absorbency level for your needs.
- ◆ Just throw your tampon away in the garbage (not in the toilet!), wash your hands, and you're done!

CAUTION!

Keep in mind that tampons can be worn at night but that it is important to remember to remove and insert a fresh tampon first thing in the morning. Never use more than one tampon at a time! Remember to remove the last tampon at the end of your period! And do not use tampons before or in between periods.

If you run into any problems or have any questions don't hesitate to reach out to us! We've been there and know what it is like. Contact us: info@yoni.care

We recommend you use a tampon with the lowest absorbency that meets your needs and that you remember to wear tampons for no longer than 8 hours! Use the following chart to select the minimum absorbency to meet your needs:

1 droplet	●○○○○○	< 6 grams
2 droplets	●●○○○○	6–9 grams
3 droplets	●●●○○○	9–12 grams
4 droplets	●●●●○○	12–15 grams
5 droplets	●●●●●○	15–18 grams
6 droplets	●●●●●●	18–21 grams

What you should know about Toxic Shock Syndrome (TSS)

TSS is a rare but serious illness that can occur among men, women and children. It is caused by toxins produced by the bacteria *Staphylococcus aureus*. This bacteria is commonly found on the skin, and inside the nose and vagina.

Approximately half of the cases of TSS are by menstruating teens and women under thirty years of age. TSS is associated with tampon use, but the disease can also occur when tampons are not used.

Symptoms

If during your period you get any of the following symptoms immediately remove your tampon and contact your doctor. The initial symptoms are similar to the flu. They usually appear very quickly and are often severe.

Symptoms include: high fever (>39°C/102°F), vomiting, diarrhea, a sunburn-like rash, dizziness and/or fainting.

Tell your doctor that you have been using tampons and that you are concerned about TSS. TSS can have fatal consequences, but can be easily treated at an early stage. It is therefore important to consult a doctor as soon as possible.

Seek medical advice before resuming the use of tampons if you have had TSS warning signs in the past. You can reduce the risk of TSS by using tampons with the lowest absorbency necessary and by changing your tampon regularly (every 4 to max. 8 hours). You can also alternate with the use of Yoni pads or panty liners instead of tampons. If you don't use tampons the risks of TSS is reduced but not entirely eliminated.

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