Bedienungsanleitung/Garantie

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Safety Precautions

Read the manual carefully before using this appliance and keep it in a safe place, together with warranty card, proof of purchase and packaging materials. If you pass the appliance on to a third party, also pass on the user manual.

- Use the appliance only for the intended, private use. This appliance is not intended for commercial use.
Do not use the appliance outdoors. Protect it from heat, direct sunlight, humidity (never submerge in any liquids) and sharp edges.
- The power adaptor is only intended for use in dry environment.
- Ensure that the output voltage, currency and polarity of the power adaptor correspond to the nameplate specifications!
- Do not touch the power adaptor with wet hands!
- Unplug the power adaptor from the appliance when cleaning or if a fault occurs.
- Make sure that the adaptor cable is not bent, crushed, run over or does not come into contact with heat sources.
- Do not operate the machine without supervision. If you leave the room you should always turn the device off. Remove the plug from the socket.
- The device and the mains lead have to be checked regularly for signs of damage. If damage is found the device must not be used.
- Do not try to repair the appliance on your own. Always contact an authorized technician. To avoid exposure to danger, always have a faulty cable replaced only by the manufacturer, by our customer service or by a qualified person and with a cable of the same type.
- Use only original spare parts.
- Do not use the appliance under a blanket.

Children and Frail Individuals

- In order to ensure your children's safety, please keep all packaging (plastic bags, boxes, polystyrene etc.) out of their reach.

WARNING!

Caution! Do not allow small children to play with the foil as there is a danger of suffocation!

- This device is not intended to be used by individuals (including children) who have restricted physical, sensory or mental abilities and/or insufficient knowledge and/or experience, unless they are supervised by an individual who is responsible for their safety or have received instructions on how to use the device.
- Children should be supervised at all times in order to ensure that they do not play with the device.

Symbols in these Instructions for Use

Important information for your safety is specially marked. It is essential to comply with these instructions in order to avoid accidents and prevent damage to the machine:

WARNING:

This warns you of dangers to your health and indicates possible injury risks.

CAUTION:

This refers to possible hazards to the machine or other objects.

NOTE:

This highlights tips and information.

Intended Use

This appliance offers a relaxing massage to adults on certain parts of the body. It cannot replace medical treatment. Do not use the massage pillow if one or more of the following warnings apply to you. If in doubt, ask your physician if the massage device is suitable for you.

WARNING:

Do not use this massage device

- after a slipped disc or pathological changes to the vertebrae.
- immediately after surgery.
- under any circumstances when you wear a pacemaker.
- if you suffer from skin diseases.
- during pregnancy.
- while sleeping.

NOTE:

- The appliance can be used on-the-move wherever a 12V socket is available, e.g. when camping.
- Do not use the device in a car.

The appliance is only intended for the purpose described in this manual. ETV Elektro-technische Vertriebsgesellschaft mbH is not responsible for any damage due to improper or negligent use.

Start-Up

- Remove all packaging materials.
- Lay the Massage Mat on a stable base (chair, armchair, bed, couch, floor).
- If necessary use the rubber loops to secure the mat.

Electrical Connection

- Connect the power adaptor to the Massage Mat with a power supply cable.
- Plug the adaptor into a properly installed 230V~ 50Hz socket. The display on the remote control is activated.

■ NOTE: 12V Adaptor

- Connect the 12V adaptor to the Massage Mat with a power supply cable.
- Plug the adaptor into a 12V socket.

On/Off

The appliance is operated via a remote control.

- POWER button: Turns the Massage function on and off.
- HEAT button: Turns the Heat function on and off.

Automatic Turn-Off function

The appliance has an automatic turn-off function after 30 minutes. It functions differently for the two operating modes (Massage / Heat).

- In pure Massage mode, the remaining time is indicated on the display by "TIME".
- If Heat is activated separately, the remaining time is not displayed.
- If the Heat function is activated later during Massage, the Heat function will only be available for the remaining massage time.
- If the Massage function is activated after the Heat function, the Heat function switches off at the same time as the massage. i.e. the total running time up to automatic turn-off is the sum of the running time of the Heat function plus 30 minutes massage time.

User Instructions**■ NOTE:**

- The massage should be comfortable and without pain. When experiencing discomfort or pain, stop operation immediately.
- Do not use the appliance
 - shortly before going to bed. The massage has a stimulating effect and so may cause restlessness.
 - for longer than 30 minutes. You can repeat the massage several times per day. Take rest breaks of at least 60 minutes between massages.
- Position yourself on the Massage Mat.
- Use the POWER button to turn on the appliance.

■ NOTE:

The appliance starts the first Massage program in the "MODE 0" setting.

- Choose between 8 different massage programs ("MODE 0" to "MODE 7"). The 6 massage zones are activated in varying sequences, intensity and rhythms. Select by pressing the MODE button on the remote control.
- Press the HEAT button on the remote control to activate the Heat function on the massage mat.

⚠ CAUTION:

People who are sensitive to heat must take particular care when using the appliance.

■ NOTE:

- You can turn the Heat function on and off at any time regardless of the Massage function.
- The activated Heat function is indicated on the display by the  symbol.

- Turn off the appliance after use.
- Unplug the mains plug (or the adaptor) from the mains outlet.

Cleaning and Storage**⚠ WARNING:**

- Before cleaning, always unplug the mains plug (or the adaptor) from the mains outlet!
- Under no circumstances should you immerse the device in water for cleaning purposes. Otherwise this might result in an electric shock or fire.

⚠ CAUTION:

- Do not use a wire brush or any abrasive items.
- Do not use any acidic or abrasive detergents.
- Use a clothes brush to remove fluff.
- Use a slightly damp cloth to remove marks. If necessary, only use a neutral cleaning agent.
- Only use the Massage Mat again when it is completely dry.
- Store in a cardboard box. Keep the appliance in a dry place.

Technical Specifications

Model: MM 5568
Net weight: 1.10kg

Power adaptor

Protection class: II
Input: AC 100-240V, 50/60Hz
Output: DC 12V, 1000mA
Polarity: 

Subject to change without prior notice!

This device complies with all current CE directives, such as electromagnetic compatibility and low voltage directive and is manufactured according to the latest safety regulations.

Meaning of the “Dustbin” Symbol

Protect our environment: do not dispose of electrical equipment in the domestic waste.

Please return any electrical equipment that you will no longer use to the collection points provided for their disposal.

This helps avoid the potential effects of incorrect disposal on the environment and human health.

This will contribute to the recycling and other forms of reutilisation of electrical and electronic equipment.

Information concerning where the equipment can be disposed of can be obtained from your local authority.

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