PPW3330

BOSCH

de Gebrauchsanleitung
en Operating instructions
fr Notice d'utilisation
it Istruzioni per l'uso
nl Gebruiksaanwijzing
da Brugsanvisning
no Bruksanvisning
sv Bruksanvisning
fi Käyttöohje
es Instrucciones de uso

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Keep the operating instructions in a safe place!

Safety instructions

This appliance is designed for the household or for similar, non-industrial and/or medical applications.

Non-industrial applications include e.g. use in employee rooms in shops, offices, agricultural and other commercial businesses, as well as use by guests in boarding houses, small hotels and similar dwellings.

⚠ Warning!

Pregnant women and persons with implanted medical devices such as pacemakers may use the weighing function only.

⚠ Risk of slipping!

Do not stand on the scales if your feet or the surface of the scales are wet – risk of slipping!

- Supervise children to prevent them from playing with the appliance.
- Do not drop anything onto the scales.
- To prevent danger, the appliance may be repaired by our customer service only.

Overview

Fig. 1

- 1 Display
- 2 Navigation buttons ▼▲
- 3 **USER** button
- 4 **SET** button
- 5 Battery compartment

Body height	100–220 cm	
Weight	kg/g (100 g graduation)	
Age	10–99 years	
Sex	Male	Female
Mode	Athlete (yes/no)	
Body Mass Index	kg/m ²	
Body fat content	as % (0.1% graduation)	
Body water content	as % (0.1% graduation)	
Muscle mass	as % (0.1% graduation)	
Bone mass	in kg/g (100 g graduation)	

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Switching on the appliance

- Insert the enclosed battery into the battery compartment **5** (Fig. 2).
- Place the scales on a level, hard and dry surface.
- To switch on the scales, tap briefly and audibly on the surface of the scales with the ball of your foot.
The vibration sensor activates the scales which require a few seconds to calibrate. When   is indicated on the display, the scales are ready.

Important!

Do not move the scales during or after calibration to ensure that the subsequent measurements are implemented exactly. If the scales are moved to a different location, repeat the process.

Use

Two functions can be selected with these scales:

- **Weighing function**
- **Analysis function**

Weighing function

Fig. 3

- To switch on the scales, tap briefly and audibly on the surface of the scales with the ball of your foot until   is indicated.
- Step onto the scales and stand still for as long as the display flickers. Distribute your weight evenly on both legs.
- The measurement can now be read off in kg/g.



The scales automatically switch off for a couple of seconds after the measuring process.

Analysis function

Important: Pregnant women and persons with implanted medical devices such as pacemakers may use the weighing function only.

To measure the percentage of body fat, body water, muscle and bone mass and to calculate your BMI, first input your personal data.

The scales have 10 memory locations in which personal data, such as size, age, sex and athlete mode, is input.

Note: If no button is pressed within 15 sec. during programming, the scales automatically switch off. If the Navigation buttons   are held down for longer, the numerical sequence is accelerated.

Inputting personal data

Fig. 4

- Switch on the scales and wait until   is indicated.
- Press the **SET** button. The setting for the personal memory location flashes.
- Using the Navigation buttons  , select the personal memory location from 1–10.
- Confirm by pressing **SET** button. The symbol for  (female) or  (male) now flashes.
- Using the Navigation buttons  , select your sex.
- Confirm by pressing **SET** button. The setting for input age flashes.
- Using the Navigation buttons  , set your age.
- Press the **SET** button. The setting for height input flashes.
- Using the Navigation buttons  , set your height.
- The symbol for Normal or Athlete mode flashes.
- Using the Navigation buttons  , select athlete yes or no.
- Confirm by pressing **SET** button.

Now all settings are displayed for approx. 2 seconds,   is indicated on the display. The scales automatically switch off after a couple of seconds.

To assign other memory locations, repeat the Inputting personal data process.

Note: Memory locations which are already occupied can be reprogrammed at any time.

Starting the analysis

Measurement of weight as well as BMI, the percentage of body fat, body water, muscle and bone mass.

Fig. 5


- Switch on the scales and wait until   is indicated.
- Using the **USER** button, select the programmed personal memory location, e.g. 2. The input personal data is displayed for 2 seconds.
- As soon as the scales display  , step on the scales and stand still for as long as the display flickers. Distribute your weight evenly on both legs.

Note: Stand on the scales with bare feet only. Initial display for approx. 3 seconds

- Weight in kg

The analysis function now requires a few seconds to perform the measurement and calculation.

Successive displays on display

- Weight in kg
- the Body Mass Index **BMI**,
- the body fat percentage **FAT**,
- the body water percentage ,
- the muscle mass  as %,
- the bone mass  in kg,

for approx. 3 seconds each.

These displays are automatically repeated 2x in succession.

Important: Stand with feet, calves and thighs apart, otherwise the measurement will be incorrect.

The scales automatically switch off a couple of seconds after the end of the measuring process.

General information

The measurements are only an indication of long-term changes. You are not replacing medical observation or consultation.

For more detailed information on the meaning of the body fat content, body water content as well as bone mass and muscle mass and the values ideal for you, please consult your GP.

Tip: To obtain comparable measurements, weigh yourself every day at the same time and under the same conditions (recommended three hours after eating or sport).

Note: The measurement results are not reliable for children under 10 years of age.

What does BMI mean?

The Body Mass Index (BMI) is a measured value for evaluating the weight of a person.

The value can be used to determine whether the person is normal weight, overweight or underweight.

The Body Mass Index (BMI) is the quotient of your weight (kg) and the square of your height (m). The unit of the BMI is kg/m².

Meaning of BMI (according to WHO)

BMI value	Meaning
under 19	underweight
19 to 25	normal weight
25 to 30	overweight
30 to 35	obesity degree 1
35 to 40	obesity degree 2
over 40	obesity degree 3

The Body Mass Index is a very approximate value only.

How does body fat measurement function?

Body fat is measured according to the BIA principle (bioelectrical impedance analysis).

The appliance transmits a weak, harmless electrical signal through the body.


The resistance (impedance), i.e. the signal which the body resists, is measured.

The resistance is affected by the composition of the human body.

Muscles and organs contain water and are therefore good conductors. Fat resists the signal and is therefore a poor conductor.

The body fat percentage is calculated from this measured resistance taking into account the data (size, age, sex) which was input into the scales.

The body fat content varies greatly from person to person. Therefore, only approximate values are indicated in the following table.

Water measured values

The body fat and water measurements are influenced by the amount of water in the body. The amount of water varies, e.g. according to the time of day, demands on the body, baths or showers, menstruation, eating and drinking. Even illnesses or medication may influence the result.

The following guide values, according to WHO, apply to persons who are normal weight:

- Women 50–55%
- Men 60–65%

What is the muscle mass?

This function indicates the calculated mass of skeletal muscles (approximate value) in your body. A human being has over 600 muscles. Muscles play an important role (movement, posture) and use energy. With increasing muscle mass (increase in activity and sport), more energy and calories are used and the body fat content is reduced.

Average skeletal muscle mass (International Commission on Radiological Protection, 1975):

- Women who weigh 58 kg: 30%
- Men who weigh 70 kg: 40%

What is bone mass?

This function indicates the calculated bone mass (approximate value for bone mineral content). The bone mass is dependent on life style. Therefore, it is important to have a healthy diet and to exercise regularly.

Please consult a specialist for the ideal values.

Average bone mineral content (Ricoet al. 1993):

- Women between 17 and 28 years of age: approx. 2.7 kg
- Men between 17 and 28 years of age: 2.7–3.5 kg

Note: For the following persons it may not be possible to calculate the bone mass.

Please use these values as an approximation.

- Elderly people or adolescents
- Women during or after the menopause
- People who are receiving hormone therapy
- Women who are breast-feeding

Athlete mode

Scientific studies show that the body composition of performance athletes differs greatly from that of amateur athletes and non-athletes. That is why these scales can be used to select between the Normal and Athlete setting. Requirement for the selection of Athlete is e.g. intensive endurance training of at least 8–10 hours per week, intensive resistance training for power athletes and a rest pulse of max. 60 beats per minute.

Cleaning the appliance

Never immerse the scales in water.

Do not use a steam cleaner.

- Wipe the scales with a damp cloth, do not use scouring agents.
- Clean the glass surface with a standard glass cleaner.

Error messages

If **Err** is indicated on the display, this may be due to various causes:

- The calibration was defective.
→ Remove and re-insert the battery. Then tap briefly and audibly on the surface of the scales with the ball of your foot.
- The weighing process was not conducted with the person standing still or for long enough. → You must stand still on the scales for at least five seconds.

- The weight to be weighed is too heavy (max. 180 kg).
- The contact between the soles of the feet and the measuring contacts is too low.
→ Place feet on the measuring contacts.

Note: Very dry skin and/or very hard skin may distort the result. The measurement results will be more accurate if the soles of the feet have good blood circulation or are damp.

However, to avoid the risk of slipping, please ensure that the soles of the feet are not wet.

- The percentage of body fat is outside the range of 5–50%
→ Repeat the measuring process
- The percentage of body water content is outside the measurement range.

If **L** or **o** is indicated, this has the following cause:

- The battery is discharged or too weak.
Change the battery (see **Fig. 2**).

Then calibrate the scales, tap briefly and audibly on the surface of the scales with the ball of your foot.

Do not use rechargeable batteries!
Use the same type of batteries only.

Disposal

This appliance has been identified in accordance with the European directive 2002/96/EG on Waste Electrical and Electronic Equipment – WEEE). The Directive paves the way for effective EU-wide withdrawal and utilization of waste appliances.

Please ask your dealer or inquire at your local authority about current means of disposal.

Dispose of old batteries at public collection centres or return to your dealer.

Guarantee

The guarantee conditions for this appliance are as defined by our representative in the country in which it is sold. Details regarding these conditions can be obtained from the dealer from whom the appliance was purchased.

The bill of sale or receipt must be produced when making any claim under the terms of this guarantee.

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Robert Bosch Hausgeräte GmbH