Technologies®  
**EMAG**  
www.emmi-dent.de

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
# EMMI-DENTAL PROFESSIONAL / EMMI-DENT 6 PROFESSIONAL

Dental Micro Cleaning with Ultrasound  
Operating Instructions – 230 V~ / 50 Hz

## Scope of Delivery

- 1 Charging Station
- 1 Handpiece
- 2 Ultrasound Attachments
- 1 Tartar Remover (Emmi-dental Professional only)
- 2 Colour Tablets (Emmi-dental Professional only – to check the difference in cleaning)
- 1 Tube of Emmi-dent Ultrasound Toothpaste

**emmi-dental**  
PROFESSIONAL

**emmi-dent 6**  
PROFESSIONAL

## Dental Micro Cleaning with Real Ultrasound It Functions so simply!

- 1** Connect the charging station with the socket.

- 2** Insert the handpiece from above into the charging station (charging time: ca. 24 hours)

- 3** Now insert the ultrasound attachment vertically from above into the handpiece and turn on the device.

- 4** Briefly wet the bristles of the ultrasound attachment under running water.

- 5** Apply only Ultrasound Toothpaste (the formation of millions of micro-bubbles only functions with this) to the bristles (approximately twice the size of a pea).

- 6** Lightly spread the toothpaste over the inside and outside of your teeth with the bristles of the ultrasound attachment.

- 7** Then hold the bristles motionless on the teeth, slightly from inside, outside and from above (chewing surface) – without applying pressure – for ca. 5 to 10 seconds. The bristles of the ultrasound attachment each cover 2-3 teeth. The actual retention time depends on the degree of impurity of the teeth.

- 8** After completion of dental cleaning (while the device is turned on), briefly rinse the bristles and handpiece under running water, turn the device off and briefly dry the handpiece (not the bristles). Ultrasound automatically disinfects the bristles.

- 9** You can simply plug the device with inserted ultrasound attachment into the charging station again after dental micro cleaning – or if you are for instance travelling, you can operate the device for ca. 10 days without recharging, if the batteries are fully charged (after ca. 24 hours of charging time).

- 10** If you clean your teeth twice a day, please replace the ultrasound attachment after 3 months.

## Intended Use

The products may only be utilised for cases of application intended in the catalogues and the technical descriptions, i.e. for cleaning teeth, and only in conjunction with external devices and components recommended or approved by EMAG AG. Any other use or additional use is not considered to be the intended use.

In the event of unintended use, EMAG AG does not assume any liability for possibly arising damages and no guarantee for flawless and correct functioning of components.

Observance of all details in the operating instructions also belongs to intended use. The  symbol on the identification plate refers to this.

The flawless and safe operation of products presupposes appropriate transport, storage, installation and assembly as well as careful operation and maintenance. The recurring tests in accordance with the German Occupational Health & Safety Regulation BGV A3 are to be heeded.

The device must be installed so that the charging station's power plug is freely accessible at any time.

## Improper Use

Dangers to life and limb as well as dangers to property can arise with improper use of the product. Improper use is prohibited for this reason. In the event of unintended use, EMAG AG does not assume any liability for possibly arising damages and no guarantee for flawless and correct functioning of components. In order to preclude an improper use of the device, the following details are to be observed:

- The ultrasonic device for dental micro cleaning may only be operated with the mains voltage which is printed on the identification plate (underside of the charging station).
- Do not use any combustible or aggressive chemicals (acids, etc.) for cleaning.
- If transport damages are detected, do not connect the charging station to the mains.
- The device may only be operated in interior rooms at room temperature.

## Operating Instructions

*IMPORTANT before initial use:*

*Before using, connect the charging station's plug with the outlet (230 V~/50 Hz), insert the handpiece into the charging station and charge for approx. 24 hours.*

While the handpiece is inserted in the charging station, the red LED display below the switch on the handpiece is illuminated in order to show that the device is flawlessly charging (display also illuminates when batteries are fully charged and the handpiece is put back on the charging station; in the process, the charging is automatically turned off). Remove the handpiece from the charging station and insert an ultrasound attachment vertically from above (please make sure that the ultrasound attachment is inserted vertically, and is also removed vertically upward in order to prevent damage to the contact pins of the handpiece; do not tilt!). Briefly wet the bristles of the ultrasound attachment under running water, and apply Ultrasound Toothpaste to the bristles (approx. twice the size of a pea).

*Please note: Use only Emmi-dent Ultrasound Toothpaste, since commercially available toothpastes do not enable a formation of microcavities (micro-bubbles).*

Turn the device on. The ultrasound is generated in the ultrasound attachment directly above the bristles and is controlled through the electronics in the handpiece. After turning the device on – without ultrasound attachment – you will feel a vibration. This means that the electronics of the handpiece are functioning flawlessly. The green LED display below the switch on the handpiece is not illuminated. With inserted ultrasound attachment – device turned on – the green LED display below the switch is illuminated and thus shows that the ultrasound generator is functioning with the proper frequency and flawlessly. If the green LED display is not illuminated, the ultrasound attachment can have a poor contact or no contact with the pins of the handpiece or can be defective. In this case, use another ultrasound attachment and contact our service department because of the defective attachment.

After turning the device on you will feel a slight vibration and hear a faint noise. This has nothing to do with the cleaning! Since one does not see, hear and feel ultrasound, a small motor has been integrated in the handpiece so that the user feels and hears that the device is turned on (as he/she is accustomed to with any other device). Lightly spread the toothpaste over the inside and outside of your teeth with the bristles (device turned on).

Then hold the bristles motionless on the teeth, slightly from inside, outside and from above (chewing surface) – without applying pressure – for ca. 5 to 10 seconds.

The bristles of the ultrasound attachment cover 2-3 teeth. The actual retention time depends on the degree of impurity of the teeth (e.g. dental plaque, nicotine, stains from wine, etc.). With acute problems, such as gingivitis or periodontitis, keep the bristles on the affected teeth with the bristles for 15 to 20 seconds. This has a faster and more efficient therapeutic effect (the duration is harmless, since this involves low-energy ultrasound). The ultrasonic cleaning of teeth usually requires ca. 4 to 6 minutes of time (perhaps a bit longer in the beginning). After completion of dental cleaning (device turned on), go a few times over the teeth in order to dissolve any detached or still adhesive impurities. Rinse the mouth as usual with water afterwards.

Briefly rinse the bristles and handpiece (while still turned on) under running water, then turn the device off and lightly dry the handpiece (not the bristles). Please make sure that the bristles are clean and free of food particles or other impurities. Since ultrasound sterilises, the bristles are automatically and thoroughly disinfected after using the device. You can simply plug the device with inserted ultrasound attachment into the charging station again after dental cleaning, or you can use the device up to ca. 10 days (with use twice a day) without recharging, if the batteries are fully charged and the daily cleaning time does not exceed a total of 10 minutes.

**Caution:**

*Absolutely make sure that the handpiece is turned off with an empty battery, since otherwise no recharging takes place! Even if the battery is not yet discharged, never insert the device into the charging station while it is turned on! Dentists recommend performing ultrasonic cleaning of teeth twice a day for ca. 3-6 minutes. If you clean your teeth twice a day, please replace the ultrasound attachment after ca. 3 months. If you take the device with you while traveling, the stored energy of the fully charged battery lasts for up to ca. 10 days, with use twice a day for ca. 2 x 5 minutes. A 24 hour bulk charge has to take place again after that.*

**Mode of Operation:**

Our ultrasonic devices are the latest and most effective method of oral hygiene and dental care (dental micro cleaning). They have been developed according to the latest scientific knowledge by German engineers, and do not function like conventional toothbrushes (dental macro cleaning), which are said to take effect through mechanical friction and abrasive substances in the toothpaste, but utilise the characteristics of modern ultrasound technologies.

Ultrasound technology has already been applied for many decades with great success in many medical fields, among other things with cerebral, ocular, pregnancy and gastric examinations, etc. The use of ultrasound for oral hygiene and for the cleaning and care of teeth is absolutely harmless for human beings, whether old or young or people with a cardiac pacemaker. It is scientifically proven that ultrasound with an output up to 50 watts is harmless to the human body. Emmi-dental functions with an output of only 0.2 watt. The effect and the safety of ultrasonic waves have been verified by dentists and implantology specialists. Many years of experience, also with crowns, inlays, implants or bridges, have not shown the slightest negative effects. The ultrasonic waves are conducted via the bristles into the Ultrasound Toothpaste medium, where they form millions of microcavities (micro-bubbles). Due to their minuteness the bubbles also penetrate into the narrowest interdental spaces, fissures, gingival pockets and clefts, e.g. with implants or braces, without any mechanical pressure and without any mechanical friction or other stress. All impurities such as tartar, food particles, bacteria etc. are thoroughly removed through implosion of micro-bubbles. At the same time, the ultrasonic waves have a positive effect on the gums, palate and tongue.

Among other things, the positive effect of ultrasonic treatment is also based on the release of the body's own energies, the promotion of lymphatic drainage and the positive influencing of the blood flow, since the ultrasonic waves penetrate into the lowest gingival layers. The prophylactic and therapeutic effect of dental cleaning and oral hygiene is based mainly on the fact that ultrasound destroys germs and bacteria, not only on the teeth and the gums, but also in the gums.

The following further purchase parts are offered for the device: *Fig. 2, page 3.*

*You can obtain the further purchase parts where you purchased your Emmi-dental device, in all pharmacies or at [www.emmi-dent.de](http://www.emmi-dent.de).*

## Ultrasound Attachments

Please note that according to statements by dentists, for hygienic reasons the ultrasound attachment is not suitable for use by several people. The same ultrasound attachment should always be used by one person. For this reason, the bristles of the ultrasound attachments are marked in colour in order to make the assignment to several family members easy. The ultrasound attachments should be replaced after ca. three months, since the bristles start to become soft from this point in time, which deteriorates the ultrasound transmission. For this purpose, remove the old ultrasound attachment vertically upward (please do not remove sideways, you can damage the contact pins of the handpiece while doing so). Firmly press the new ultrasound attachment vertically from above into the handpiece, with bristles facing forward in the direction of the start button (caution: there is a snap lock; please carefully press open, do not pinch your skin while doing so).

## Ultrasound Toothpaste

The erroneous opinion is often expressed that ultrasound requires water as a transmission medium (as is frequently known from ultrasonic baths). It is correct that the suitable transmission medium depends on the ultrasonic frequency and ultrasonic output. Examples are generally known from other medical applications such as cerebral examinations, ocular examinations, gastric examinations with ultrasound, etc. The suitable transmission media, ointments and gels, are of different consistency here. A special toothpaste (Emmi-dent) is used for dental micro cleaning, without which the ultrasound cannot take effect (millions of micro-bubbles only form in this medium). The ultrasound cannot form any micro-bubbles in conventional toothpastes. One reason for this is that these toothpastes contain a relatively high percentage of "abrasive substances" (a measurement for this is the so-called RDA value, which lies at 70 – 130 with conventional toothpastes; in comparison, the RDA value with dental micro cleaning is nil). Therefore, please note that Emmi-dent Ultrasound Toothpaste must be used for dental micro cleaning.

## Tartar Remover

The Tartar Remover functions mechanically without ultrasound. It can, but does not have to be, inserted into the handpiece. It removes strong tartar by carefully abrading the dental surfaces with the dark tip (please do not rub too strongly so as not to damage the dental enamel).

## Battery Care and Disposal

The batteries are in the device's waterproof handpiece. After use, the handpiece can be placed back in the charging station again, or if it is fully charged, it can be used for up to ca. 10 days without recharging. If the total cleaning time lasts longer than 10 minutes (e.g. with several users), the batteries can be recharged in an appropriately shorter period of time. Absolutely make sure that the batteries are never completely discharged. Once the batteries have reached a critical discharge value, the green display is not illuminated when the device is turned on. The display is illuminated in red and starts to blink after 5 seconds. If this is the case, the device absolutely has to be recharged, since otherwise the batteries can be damaged.

The device contains special environmentally-friendly NiMH batteries which may not be disposed of with household waste. In Germany, the collection containers of the GRS Batteries Foundation (common collection system) and other collection points (e.g. all points of sale for similar batteries) are available for disposal).

## GUARANTEE

There is a guarantee of 24 months – from date of purchase – for this product. Defects which arise during this guarantee period in the form of material and/or manufacturing defects will be remedied at no charge, either by means of repair or through replacement delivery. Insofar as legally permissible, any other claims for compensation are excluded. The guarantee is inapplicable in the event of intervention by third parties or disassembly by third parties without our prior consent. The guarantee is also inapplicable in the event of arbitrary damage or incorrect handling. The guarantee period is not extended with provision of a guarantee. We reserve the right to invoice the costs incurred for unjustified complaints such as installation or operating errors. For making a guarantee claim it is necessary to save the proof of purchase with date of purchase.

For fast handling of your guarantee claim cases, please send these – irrespective of specialised dealer – post-paid directly to EMAG AG. Unfortunately, unpaid shipments cannot be accepted due to logistical reasons.

Thank you for your cooperation and understanding.

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Date of Sale

Vendor's Stamp

Copyright EMAG AG. All rights reserved. Manufactured in FRG / EU. The contents of these operating instructions are subject to copyright and may neither be entirely or partially copied. Changes, abridgements, enhancements and additions, any publication or translation by third parties requires prior, written consent by EMAG AG.

# EMMI-DENTAL PROFESSIONAL / EMMI-DENT 6 PROFESSIONAL

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© Copyright EMAG AG. Alle Rechte vorbehalten. Hergestellt in der BRD/EU. Der Inhalt dieser Bedienungsanleitung unterliegt dem Urheberrecht und sie darf weder ganz noch teilweise kopiert werden. Veränderungen, Kürzungen, Erweiterungen und Ergänzungen jede Veröffentlichung oder Übersetzung durch Dritte bedarf der vorherigen schriftlichen Einwilligung durch EMAG AG.
Copyright EMAG AG. Alle rechten voorbehouden. Geproduceerd in de BRD / EU. De inhoud van deze gebruiksaanwijzing is onderworpen aan het auteursrecht en mag noch geheel, noch gedeeltelijk worden gekopieerd. Voor wijzigingen, inkortingen, uitbreidingen en aanvullingen, iedere publicatie of vertaling door derden is de voorafgaande, schriftelijke toestemming door EMAG AG vereist.
Copyright EMAG AG.. Tutti i diritti riservati. Prodotto nella RFT/UE. Il contenuto delle presenti Istruzioni per l'uso è protetto da diritto d'autore e non può essere copiato integralmente o parzialmente. Eventuali modifiche, riduzioni, ampliamenti e integrazioni, qualsiasi pubblicazione o traduzione da parte di terzi sono soggette all'approvazione scritta da parte di EMAG AG.