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


# ENGLISH

Instruction manual Inventum bathroom scale, PW821BG.

Please read these instructions carefully and completely before operating this appliance and retain for future use.

Maximum capacity: 150 kg/24 St/ 330 lb

Minimum capacity: 10 kg/1 St 6/ 20 lb

Division: per 100 gram/ ¼ St/ 0,2 lb

- A Electrodes
- B Platform
- C Down-key
- D SET-key (confirm)

- E Up-key
- F Display
- G Convert-switch (Kg/Lb/St)

## EXPLANATION OF SYMBOLS

Legislation requires that all electrical and electronic equipment must be collected for reuse and recycling. Electrical and electronic equipment marked with the symbol indicating separate collection of such equipment must be returned to a municipal waste collection point.

## SAFETY INSTRUCTIONS

### • Please read the instructions carefully and completely before operating this appliance and retain for future use.

- The appliance is intended solely for domestic use within the home.
- Always use batteries which corresponds with the voltage printed on the bottom of the appliance.
- **Attention:** Do not use during pregnancy. Inaccurate measurements may be caused by the amniotic fluid.
- **Attention: The scale may not be used by persons with medical implants (e.g. pacemakers). Otherwise their function may be impaired.**
- If possible, always weigh yourself at the same time of day (preferably in the morning), after using the toilet, on an empty stomach and unclothed, in order to achieve comparable results.
- Stand on the scale without moving and with equal weight distribution on both feet!
- Minor deviations are normal. When you drink 1 glass of water for example a difference of 300 grams can occur.
- Caution! Do not stand on the scale with wet feet and do not step on the scale when its surface is wet – danger of slipping!
- Protect the scale against hard knocks, moisture, dust, chemicals, major temperature fluctuations and heat sources which are too close (stove, heating radiators).
- When you want to store the scale for a longer period, we advise to keep it in the original carton.
- Always ensure that the scale is placed on a flat, stable surface and only use it indoors.

- Repairs may only be carried out by Inventum customer service or authorized dealers. Before submitting any complaint, first check the batteries and replace them if necessary.
- If you decide never to use the appliance again, we advise you to dispose of it in special containers or at a municipal waste collection point.
- **Never immerse the scale in water or other liquids.**
- The device is only for private use and is not intended for medical or commercial purposes.
- The scale can accept a maximum load of 150 kg (330 lb, 24St). The weight measurement is displayed in increments of 100 g (0.2lb, 1/4St). The measuring results for the body fat, body water and muscle percentages are displayed in 0.1% increments.
- On delivery the scale is set to the units "cm" and "kg". A switch located on the back of the scale enables you to change to "inch", "pound" and "stones" (lb, St).
- If the appliance is used for other purposes than intended, or is it not handled in accordance with the instruction manual, the full responsibility for any consequences will rest with the user. Any damages to the product or other things are not covered by the warranty.

## INSTRUCTIONS FOR USE

### **Interesting facts**

#### **The measuring principle of the diagnostic scale**

This scale determines the body weight, bodyfat, body water and the percentage of muscle according to the principle of Bioelectric Impedance Analysis (BIA). With this process, an imperceptible, completely safe and harmless current enables the determination of body tissue proportions within seconds. This measurement of electrical resistance (impedance) and the inclusion of constants and/or individual values (age, height, sex, activity levels) can be used to determine the percentage of fat, water and muscle in the body. Muscle tissue and water have good electrical conductivity, as the fat cells hardly conduct the current due to their very high resistance.

#### **Weighing correctly**

If possible, always weigh yourself at the same time of day (ideally in the morning), after going to the toilet, on an empty stomach and without clothes to achieve comparable results. Only the long-term trend counts.

Stand on the scale without moving and with equal weight distribution on both feet!

Minor deviations are normal. When you drink 1 glass of water for example a difference of 300 grams can occur.

Do not stand on the scale with wet feet and do not step on the scale when its surface is wet – danger of slipping! Place the scale on a firm, flat surface; a firm floor covering is a basic requirement for correct weight measurement.


### **Before operation**

Before you use the appliance for the first time, please do as follows: carefully unpack the appliance and remove all the packaging material. Keep the material (plastic bags and cardboard) out of reach of children. Check the appliance after unpacking for any damage, possible from transportation. If present, pull the battery insulating strip off the battery compartment cover or remove the battery's protective film and insert the battery according to the polarity.

### **Batteries**

Your scale is equipped with a replace battery indicator. When the scale is operated with a battery which is too weak, "Lo" appears in the display field and the scale automatically switches off. The battery must be replaced in this case (2 x 3V lithium battery, CR2032).

If the scale fails to operate, remove the battery completely and insert it again. The used, completely drained standard and rechargeable batteries must be disposed of in specially marked collection containers, at toxic waste collection points or electrical product retailers. You are legally obligated to dispose of the batteries.

**Note:** Batteries containing pollutants are marked with the following symbols: Pb = Battery contains lead, Cd = Battery contains cadmium, Hg = Battery contains mercury.

### **Measuring weight only**

- Quickly and firmly tap the standing surface of your scale with your foot. The entire display (Fig. 1) appears up to "0.0" (Fig. 2) as a self-test (tap-on technology).
- Now the scale is ready to measure your weight.
- Stand on the scale without moving about and distribute your weight evenly on both legs.
- The scale immediately begins to measure your weight.
- Your measured weight is displayed shortly after this (Fig. 3). The scale switches off a few seconds after you step down from it.

Fig. 1

Fig. 2

## **Measuring weight, body fat, body water and muscle percentage**

To calculate your body fat and other values, you must enter your personal settings. The scale has 10 memory presets for basic data with which the users, e.g. each member of the family, can store and call up their personal settings at any time.

- First you personal settings must be stored.
- Briefly step onto the scale first and wait until the display "0.0" (Fig. 2) appears. Then press "SET". Now the first memory preset flashes in the display (Fig. 4). Select the desired preset by pressing the buttons "up" ► and "down" ◀ and confirm the desired preset with the "SET" button.
- Now height in cm (inch) flashes in the display (Fig. 5). Set the desired value by briefly pressing or holding the buttons "up" ► or "down" ◀ and confirm it by pressing "SET".
- "Age" in years then flashes in the display (Fig. 6). Set the desired value by briefly pressing or holding the buttons "up" ► or "down" ◀ and confirm it by pressing "SET".
- The sex setting "Woman" or "Man" appears on the display (Fig. 7). Select your sex by pressing the buttons "up" ► or "down" ◀ and confirm by pressing "SET".
- Now the setting for your individual level of activity or fitness appears on the display (Fig. 8). Select the desired level of activity by pressing the buttons (Fig. 8). Select the desired level of activity by pressing the buttons "up" ► or "down" ◀ and confirm the desired level of activity with the "SET" button.

The values you've just set are displayed again consecutively. After which the scale switches off automatically.

**Note:** To set a memory position for an other user follow the same steps as written previously and choose the position you wish to overwrite.

Fig. 3

Fig. 4

Fig. 5

Fig. 6

Fig. 7

Fig. 8

### **Explanation of the 5 levels of activity**

When selecting the level of activity the average and long-term counts.

- **A 1:** No physical activity.

- **A 2:** Little physical activity.

Little and light physical exercises (e.g. walking, light gardening and gymnastics). This mode should also be selected before the start of a training or dietary program. After approx. 6-10 weeks you can switch to the next level of activity.

- **A 3:** Moderate physical activity.

30 minutes of physical exercise 2 to 4 times a week (e.g. jogging, cycling, tennis, etc.). This mode should be selected as soon as you notice a general improvement in your well-being/fitness in the course of a training/dietary program. After 8-12 weeks of continuous training/diet, you can change to the next level of activity.

- **A 4:** Substantial physical activity.

30 minutes of physical exercise 4 to 6 times a week.

- **A 5:** Intensive physical activity.

Daily intensive training or physical work (e.g. intensive running, construction work, etc.). For this mode you should have an overall high level of mobility, endurance and power and maintain this level over a longer period of time.

### **To carry out a measurement**

After you've entered your personal setting, weight, body fat and other values can be calculated.

- Briefly step onto the scale first and wait until the display "0.0" appears.
- Select the memory preset where your basic personal data are stored by pressing the buttons "up" ► and "down" ◀ repeatedly. The data are displayed consecutively until "0.0" (Fig. 2) appears.
- Step onto the scale barefoot and make sure you're standing on both electrodes. First your body weight is determined and displayed.
- Remain standing on the scale; now the body fat and body water analysis is performed. This can take a few seconds.

**Important:** Your feet, legs, calves and thighs must not touch each other, as otherwise the measurement cannot be performed correctly.

The following data are displayed:

- Body fat percentage in % (Fig. 9)
- Water percentage in % (Fig. 10)
- Muscle percentage in % (Fig. 11)
- Now all measured values are displayed consecutively and the scale switches off.

Fig. 9

Fig. 10

Fig. 11

## **Body fat guide**

The following body fat levels provide you with a guideline (for further information, please consult your doctor!).

### **Women**

### **Men**

Age	Very good	Good	Average	Poor	Age	Very good	Good	Average	Poor
< 19	< 17%	17-22%	22,1-27%	> 27,1%	< 19	< 12%	12-17%	17,1-22%	> 22,1%
20 - 29	< 18%	18-23%	23,1-28%	> 28,1%	20-29	< 13%	13-18%	18,1-23%	> 23,1%
30 - 39	< 19%	19-24%	24,1-29%	> 29,1%	30-39	< 14%	14-19%	19,1-24%	> 24,1%
40 - 49	< 20%	20-25%	25,1-30%	> 30,1%	40-49	< 15%	15-20%	20,1-25%	> 25,1%
> 50	< 21%	21-26%	26,1-31%	> 31,1%	> 50	< 16%	16-21%	21,1-26%	> 26,1%

based on: "Principles + Labs for Physical Fitness and Wellness. 1st edition Copyright 1999"

A lower level is often found in athletes. Depending on the type of sports, the intensity of training and the person's physical constitution, levels can be achieved that are even lower than the specified guidelines.

## **Body water percentage**

According to the World Health Organization (WHO 2001), the percentage of body water should lie within the following ranges!

Women: 50 - 55%

Men: 60 - 65%

Children: 65 - 75%

In persons with a high percentage of body fat, the percentage of body water very frequently lies below the specified guidelines. As the percentage of muscle differs from individual to individual, there are no generally valid guidelines.

## **Muscle percentage**

The percentage of muscle differs from individual to individual, they should lie within the following ranges!

### **Women**

### **Men**

Age	little	normal	many	Age	little	normal	many
< 19	< 35%	35 - 41%	> 41%	< 19	< 43%	43 - 56%	> 56%
20-29	< 34%	34 - 39%	> 39%	20-29	< 42%	42 - 54%	> 54%
30-39	< 33%	33 - 38%	> 38%	30-39	< 41%	41 - 52%	> 52%
40-49	< 31%	31 - 36%	> 36%	40-49	< 40%	40 - 50%	> 50%
> 50	< 29%	29 - 34%	> 34%	> 50	< 39%	39 - 48%	> 48%

## **Temporary relation of the results**

Only the long-term trend counts. Short-term weight fluctuations within a few days are usually only due to a loss of liquid.

The interpretation of the results depends on the changes in the total weight, percentage of fat, percentage of water, percentage of muscle and on the length of time during which these changes take place. Rapid changes over several days should be distinguished from mid-term changes (over several weeks) and long-term changes (months).

As a basic rule, short-term weight changes almost exclusively represent changes in water content, while mid-term and long-term changes can also concern the percentage of fat and muscle.

- If you experience a short-term drop in weight, however the percentage of body fat increases or remains the same, you have only lost water, e.g. after a training session, going to the sauna or following a diet limited to rapid weight loss.
- If you experience a mid-term weight increase, while the percentage of body fat drops or remains the same, you might have built up your muscle mass or your body may have retained more water (e.g. as the result of taking medication or a weaker cardiac function in elderly people).
- If your weight and percentage of body fat drop simultaneously, your diet is working and you are losing fat mass. Ideally, you can support your diet with physical exercise, fitness or power training. This can result in a mid-term increase in the percentage of muscle in your body.

The percentage body fat and body water and the percentage of muscle can not be added (muscle tissue also contains body water).

## **Tips for using the scale**

Important when measuring body fat/body water/muscle percentage:

- The measurement may only be carried out while barefoot and it is helpful if the soles of your feet are slightly damp. Completely dry soles can result in unsatisfactory results, as they have insufficient conductivity.
- Stand still during the measurement.
- Wait several hours (6-8) after unusually strenuous activity.
- Wait approx. 15 minutes after getting out of bed so that the water in your body can be distributed.

The measurement is not reliable for:

- Children under approx. 10 years of age;
- Professional sporters and bodybuilders;
- Pregnant women;
- Persons with fever, undergoing dialysis, with symptoms of edema or osteoporosis;
- Persons taking cardiovascular medication. Persons taking vasodilating or vasoconstricting medications;
- Persons with substantial anatomical deviations in the legs relative to their total height (leg length considerably shorter or longer than usual).

## **Incorrect measurement**

If the scale detects an error during weighing, «Err» appears in the display.

Possible causes of errors:	Remedy:
The scale was not activated before stepping on. When you stand on the scale before "0.0" appeared in the display, it will not operate correctly.	Activate the scale the right way and repeat the measurement.
The maximum load-bearing capacity of 150 kg. (330 lbs) was exceeded.	Only weigh the maximum permissible weight.
The electrical resistance between the electrodes and the soles of your feed is too high (e.g. with heavily callused skin).	Repeat weighing barefoot. Slightly moisten the soles of your feet if necessary. Remove the calluses on the soles of your feet if necessary.
You were not standing still on the scale.	Stand still during the measurement.
Your body fat lies outside the measurable range (less than 5% or greater than 50%).	Repeat weighing barefoot and slightly moisten the soles of your feet if necessary.
Your body water lies outside the measurable range (less than 43% or greater than 75%).	Repeat weighing barefoot and slightly moisten the soles of your feet if necessary.
The percentage of muscle lies outside the measurable range (depending on age and sex).	Repeat weighing barefoot and slightly moisten the soles of your feet if necessary.

And remember: Physical exercise is the basis for a healthy body.

## **Convert from Kg (kilogram)/ St (Stone)/ Lb (Pounds)**

This scale can display measurements in kilograms, Stone and pounds. You can set this with the switch on the bottom of the scale.  
Kg (kilograms), St (Stone) or Lb (pounds).

## C L E A N I N G & M A I N T E N A N C E

Clean the scale occasionally with a damp cloth. Do not use aggressive cleaning agents and never submerge the scale. **Never submerge the scale in water or other liquids.**

## G U A R A N T E E & S E R V I C E

- This product is guaranteed for a period of 2 years from the date of purchase against mechanical and electrical defect.
- Guarantee can only be given if proof of purchase (e.g. till receipt) can be shown. Therefore keep it in a safe place. The guarantee period starts at date of purchase.
- Guarantee is not valid in case of:
  - normal wear;
  - damage through misuse use or when the instructions were not followed correctly;
  - lack of maintenance;
  - not following the operation- and safety instructions;
  - incompetent montage or repairs made by yourself or others;
  - using accessories and spare parts which are not suitable for the appliance;for assessment by our technical department.
- The guarantee is not valid for:
  - batteries;for assessment by our technical department.

Return the appliance, including proof of purchase and description of the complaint, stamped to:

Inventum Holland BV  
Generatorstraat 17  
3903 LH VEENENDAAL  
The Netherlands

or to:  
new E.T.P. bvba  
Passtraat 121  
9100 Sint Niklaas  
BELGIUM

- Pack it carefully (preferably in the original carton). Ensure unit is clean.
- Enclose your name and address and quote model number PW821G on all correspondence.

When the guarantee is no longer valid or when the period of guarantee has passed, we will charge any costs that flow out of an assessment and/or repairing.

On request we will supply a statement before repairing. Without requesting a statement, we assume that you, by sending the appliance, automatically order to repair the appliance.

A guarantee period of two months is valid for repairs done by our technical department and the used spare parts, counting from the repairdate. When a malfunction occurs within this period which directly flows out of the repairs and the used spare parts, the consumer has the right of repairs without any costs.

We reserve the right to change this product, packaging and documentation without notice.

Inventum Holland BV  
Generatorstraat 17, 3903 LH , Veenendaal, Nederland  
Tel.: 0318-542222, Fax.: 0318-542022, Mail: [helpdesk@inventum.eu](mailto:helpdesk@inventum.eu)

Inventum Belgium BVBA  
Xavier de Cocklaan 66 bus 8, 9830, Deurle, België  
Tel: 09 220 30 20, Fax.: 09 220 30 22, Mail: [info@inventumbelgium.be](mailto:info@inventumbelgium.be)

Website: [www.inventum.eu](http://www.inventum.eu)

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Op basis van: "Principles + Labs for Physical Fitness and Wellness. 1e editie Copyright 1999"
basiert auf: "Principles + Labs for Physical Fitness and Wellness. 1st edition by Copyright 1999"
basé sur: "Principles + Labs for Physical Fitness and Wellness. 1st edition by Copyright 1999"