INVENTUM

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


ENGLISH

Instruction manual Inventum bathroom scale, PW890BG.

Please read these instructions carefully and completely before operating this appliance and retain for future use.

Maximum capacity: 150 kg/24 St/ 330 lb
Minimum capacity: 10 kg/1 St 6/ 20 lb
Division: per 100 grams/ 1/4 St/ 0,2 lb

- | | | | |
|---|-------------------|---|---------------------------|
| A | Electrodes | E | USER-key |
| B | Glass platform | F | Up-key |
| C | Down-key | G | Display |
| D | SET-key (confirm) | H | Convert switch (Kg/St/Lb) |

EXPLANATION OF SYMBOLS

Legislation requires that all electrical and electronic equipment must be collected for reuse and recycling. Electrical and electronic equipment marked with the symbol indicating separate collection of such equipment must be returned to a municipal waste collection point.

SAFETY INSTRUCTIONS

- **Please read the instructions carefully and completely before operating this appliance and retain for future use.**
- The appliance is intended solely for domestic use within the home.
- Always use batteries which corresponds with the voltage printed on the bottom of the appliance.
- **Attention:** Do not use during pregnancy. Inaccurate measurements may be caused by the amniotic fluid.
- **Attention: The scale may not be used by persons with medical implants (e.g. pacemakers). Otherwise their function may be impaired.**
- If possible, always weigh yourself at the same time of day (preferably in the morning), after using the toilet, on an empty stomach and unclothed, in order to achieve comparable results.
- Stand on the scale without moving and with equal weight distribution on both feet!
- Minor deviations are normal. When you drink 1 glass of water for example a difference of 300 grams can occur.
- Caution! Do not stand on the scale with wet feet and do not step on the scale when its surface is wet – danger of slipping!
- Protect the scale against hard knocks, moisture, dust, chemicals, major temperature fluctuations and heat sources which are too close (stove, heating radiators).
- When you want to store the scale for a longer period, we advise to keep it in the original carton.
- Always ensure that the scale is placed on a flat, stable surface and only use it indoors.

- Repairs may only be carried out by Inventum customer service or authorized dealers. Before submitting any complaint, first check the batteries and replace them if necessary.
- If you decide never to use the appliance again, we advise you to dispose of it in special containers or at a municipal waste collection point.
- **Never immerse the scale in water or other liquids.**
- The device is only for private use and is not intended for medical or commercial purposes.
- The scale can accept a maximum load of 150 kg (330 lb, 24St). The weight measurement is displayed in increments of 100 g (0.2lb, ¼St). The measuring results for the body fat, body water and muscle percentages are displayed in 0.1% increments.
- On delivery the scale is set to the units "cm" and "kg". A switch located on the back of the scale enables you to change to "inch", "pound" and "stones" (lb, St).
- If the appliance is used for other purposes than intended, or is it not handled in accordance with the instruction manual, the full responsibility for any consequences will rest with the user. Any damages to the product or other things are not covered by the warranty.

INSTRUCTIONS FOR USE

Interesting facts

The measuring principle of the diagnostic scale

This scale determines the body weight, bodyfat, body water and the percentage of muscle according to the principle of Bioelectric Impedance Analysis (BIA). With this process, an imperceptible, completely safe and harmless current enables the determination of body tissue proportions within seconds. This measurement of electrical resistance (impedance) and the inclusion of constants and/or individual values (age, height, sex, activity levels) can be used to determine the percentage of fat, water and muscle in the body. Muscle tissue and water have good electrical conductivity, as the fat cells hardly conduct the current due to their very high resistance.

Weighing correctly

If possible, always weighs yourself at the same time of day (ideally in the morning), after going to the toilet, on an empty stomach and without clothes to achieve comparable results. Only the long-term trend counts.

Stand on the scale without moving and with equal weight distribution on both feet!

Minor deviations are normal. When you drink 1 glass of water for example a difference of 300 grams can occur.

Do not stand on the scale with wet feet and do not step on the scale when its surface is wet – danger of slipping! Place the scale on a firm, flat surface; a firm floor covering is a basic requirement for correct weight measurement.

Before operation

Before you use the appliance for the first time, please do as follows: carefully unpack the appliance and remove all the packaging material. Keep the material (plastic bags and cardboard) out of reach of children. Check the appliance after unpacking for any damage, possible from transportation. If present, pull the battery insulating strip off the battery compartment cover or remove the battery's protective film and insert the battery according to the polarity.

Batteries

If provided, pull out the battery insulating strip at the battery compartment or remove the protective film from the battery and insert the battery with the correct polarity. If the scale fails to operate, remove the battery completely and reinsert it. Your scale is provided with a "Change battery" indicator. If the scale is operating with a battery which is too weak, "Lo" appears on the display and the scale will switch off automatically. In this case, the batteries must be exchanged (3 x 1,5 V, AAA, Mikro). The used, completely drained standard and rechargeable batteries must be disposed of in specially marked collection containers, at toxic waste collection points or electrical product retailers. You are legally obligated to dispose of the batteries.

Note: Batteries containing pollutants are marked with the following symbols: Pb = Battery contains lead, Cd = Battery contains cadmium, Hg = Battery contains mercury.

Measuring weight only

- Press quickly and forcefully with your foot on the platform of your scale to activate it. As a self-test, the complete display will appear up to "0.0".
- The scale is now ready to measure your weight.
- Step on the scale and stand still with your weight distributed evenly on both legs.
- The scale will start measuring your weight immediately and a short time later will display your weight.
- When you step off the scale, the scale will switch off after a few seconds.

Measuring weight, body fat, body water and other values

To measure your body fat percentage and other physical data, you must enter your personal user parameters.

The scale has 10 memory positions in which you and other members of your family can save and recall personal settings.

- Switch on the scale by briefly stepping on the platform. Wait until the display shows "0.0". Then press "SET". The first memory position will appear flashing on the display. Now you can enter the following settings:

Memory position	1-10
Height	100 to 220 cm (3'-03" to 7'-03")
Age	10 to 100 years
Gender	man (♂), woman (♀)
Activitylevel	1 to 5

You can enter the relevant settings with short or long pressure on the up button ▲ or down button ▼. Confirm the settings in each case by pressing "SET". The values stored in this way are then shown again consecutively, after which the scale switches off automatically.

Note: To set a memory position for an other user follow the same steps as written previously and choose the position you wish to overwrite.

Explanation of the 5 levels of activity

When selecting the level of activity the average and long-term counts.

- A 1:** No physical activity.
- A 2:** Little physical activity.

Little and light physical exercises (e.g. walking, light gardening and gymnastics). This mode should also be selected before the start of a training or dietary program. After approx. 6-10 weeks you can switch to the next level of activity.

- A 3:** Moderate physical activity.

30 minutes of physical exercise 2 to 4 times a week (e.g. jogging, cycling, tennis, etc.). This mode should be selected as soon as you notice a general improvement in your well-being/fitness in the course of a training/dietary program. After 8-12 weeks of continuous training/diet, you can change to the next level of activity.

- A 4:** Substantial physical activity.

30 minutes of physical exercise 4 to 6 times a week.

- A 5:** Intensive physical activity.

Daily intensive training or physical work (e.g. intensive running, construction work, etc.). For this mode you should have an overall high level of mobility, endurance and power and maintain this level over a longer period of time.

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After all parameters have been entered, you can now measure your weight, body fat and the other data.

- Switch on the scale by briefly stepping on the platform. Wait until the display shows "0.0".
- Select the memory position in which your personal parameters are stored by pressing the "User" button several times. These parameters will then appear in turn until the "0.0" display appears.
- Step on the scale with bare feet and stand still on the electrodes. **Important:** There may be no contact between feet, legs, calves and thighs. Otherwise the measurement will not be performed correctly.

The following display then appears:

1. Weight, in kg.	KG
2. Body fat percentage, in %	BF
3. Body water percentage in %	≡
4. Muscle percentage in %	⚖️
5. Basal metabolic rate in Kcal	♂️

The display then changes from muscle percentage and basal metabolic rate to:

6. Bone mass in kg	⚖️
7. Active metabolic rate in Kcal	♂️

Pressing any of the buttons will show the results more rapidly.

Up to 5 results are automatically saved in the data memory. Incorrect measurements indicated with "Err" are not saved.

Measured data recall

The last 5 measurements can be recalled.

- Switch on the scale by briefly stepping on the platform.
- Select the memory position in which your personal parameters are stored by pressing the "User" button several times. These then appear in turn until the "0.0" display appears.
- As soon as "0.0" appears, you can recall the last data measured with the up button ▲.
- Press the up button ▲ repeatedly to display the next older data. After a maximum of 5 measurements, the scale returns to the "0.0" display and is ready for measurement.

If user settings are changed (e.g. age), the measured data memory is deleted.

Evaluation of results

Trend display (Trackersystem)

This scale has a trend display, in which the value currently measured is compared with the averages of the last 5 measurements.

- Trend display up ▲ : present value is higher than the average of the last 5 measurements.
- Trend display down ▼ : present value is lower than the average of the last 5 measurements.
- Trend display up and down ▲ ▼ : present value is the same as the average of the last 5 measurements.

Note that a trend is shown only for the body weight and the body fat, body water and muscle percentages and can be shown only when you have stored at least 5 body fat percentages etc.

Interpretation

The weight values and the body fat, body water and muscle percentages are interpreted with a bar in the evaluation graphic. When the values are low, the bar is shown to the left (-). When the bar is in the middle area, the values are normal. While in the right-hand area (+) the values are relatively high.

Weight:

Your weight is compared with the BMI (Body Mass Index). It is calculated using the following formula:

$$\text{BMI} = \text{Weight in kg} / (\text{height in cm})^2$$

Underweight	Normal weight	Overweight
< 20	20 - 25	> 25

Note that with very muscular physiques (body builders), the BMI interpretation will show overweight. The reason for this is that no account is taken of well above average muscle mass in the BMI formula.

Body fat percentage

The following body fat percentages are for guidance (contact your physician for further information).

Man

Age	very good	good	average	poor
10-14	< 11%	11-16%	16,1-23%	> 23,1%
15-19	< 12%	12-17%	17,1-22%	> 22,1%
20-29	< 13%	13-18%	18,1-23%	> 23,1%
30-39	< 14%	14-19%	19,1-24%	> 24,1%
40-49	< 15%	15-20%	20,1-25%	> 25,1%
50-59	< 16%	16-21%	21,1-26%	> 26,1%
60-69	< 17%	17-22%	22,1-27%	> 27,1%
70-100	< 18%	18-23%	23,1-28%	> 28,1%

Woman

Age	very good	good	average	poor
10-14	< 16%	16-21%	21,1-26%	> 26,1%
15-19	< 17%	17-22%	22,1-27%	> 27,1%
20-29	< 18%	18-23%	23,1-28%	> 28,1%
30-39	< 19%	19-24%	24,1-29%	> 29,1%
40-49	< 20%	20-25%	25,1-30%	> 30,1%
50-59	< 21%	21-26%	26,1-31%	> 31,1%
60-69	< 22%	22-27%	27,1-32%	> 32,1%
70-100	< 23%	23-28%	28,1-33%	> 33,1%

A lower value is often found in athletes. Depending on the type of sports, training intensity and physical constitution, values may result which are below the recommended values stated. It should, however, be noted that there could be a danger to health in the case of extremely low values.

Body water

The body water percentage is normally within the following ranges:

Man

Age	poor	good	very good
10-100	< 50%	50-65%	> 65%

Woman

Age	poor	good	very good
10-100	< 45%	45-60%	> 60%

Body fat contains relatively little water. Therefore persons with a high body fat percentage have body water percentages below the recommended values. With endurance athletes, however, the recommended values could be exceeded due to low fat percentages and high muscle percentage.

Body water measurement with this scale is not suitable for drawing medical conclusions, for example concerning age-related water retention. If necessary ask your physician. Basically, a high body water percentage should be the aim.

Muscle percentage

The muscle percentage is normally within the following ranges:

Man

Age	low	normal	high
10-14	< 44%	44 - 57%	> 57%
15-19	< 43%	43 - 56%	> 56%
20-29	< 42%	42 - 54%	> 54%
30-39	< 41%	41 - 52%	> 52%
40-49	< 40%	40 - 50%	> 50%
50-59	< 39%	39 - 48%	> 48%
60-69	< 38%	38 - 47%	> 47%
70-100	< 37%	37 - 46%	> 46%

Woman

Age	low	normal	high
10-14	< 36%	36 - 43%	> 43%
15-19	< 35%	35 - 41%	> 41%
20-29	< 34%	34 - 39%	> 39%
30-39	< 33%	33 - 38%	> 38%
40-49	< 31%	31 - 36%	> 36%
50-59	< 29%	29 - 34%	> 34%
60-69	< 28%	28 - 33%	> 33%
70-100	< 27%	27 - 32%	> 32%

Bone mass:

Like the rest of our body, our bones are subject to the natural development, degeneration and ageing processes. Bone mass increases rapidly in childhood and reaches its maximum between 30 and 40 years of age. Bone mass reduces slightly with increasing age. You can reduce this degeneration somewhat with healthy nutrition (particularly calcium and vitamin D) and regular exercise. With appropriate muscle building, you can also strengthen your bone structure.

Note that this scale will not show you the calcium content of your bones, but will measure the weight of all bone constituents (organic substances, inorganic substances and water).

Attention: Please do not confuse bone mass with bone density. Bone density can be determined only by means of a medical examination (e.g. computer tomography, ultrasound). It is therefore not possible to draw conclusions concerning changes to the bones and bone hardness (e.g. osteoporosis) using this scale.

Little influence can be exerted on bone mass, but it will vary slightly within the influencing factors (weight, height, age, gender).

BMR = Basal Metabolic Rate /

The basal metabolic rate (BMR) is the amount of energy required by the body at complete rest to maintain its basic functions (e.g. while lying in bed for 24 hours). This value largely depends on weight, height and age. It is displayed on the diagnostic scale in kcal/day units using the scientifically recognized Harris-Benedict formula.

Your body requires this amount of energy in any case and it must be reintroduced into your body in the form of nutrition. If you take on less energy over the longer term, this can be harmful to your health.

AMR = Active Metabolic Rate

The energy consumption of a human being rises with increasing physical activity and is measured on the diagnostic scale in relation to the degree of activity entered (1-5). To maintain your existing weight, the amount of energy used must be reintroduced into the body in the form of food and drink. If less energy is introduced than is used over a longer period of time, your body will obtain the difference largely from the amount of fat stored and your weight will reduce. If, on the other hand, over a longer period of time more energy is introduced than the total active metabolic rate (AMR) calculated, your body will be unable to burn off the excess energy, and the excess will be stored in the body as fat and your weight will increase.

Results in relation to time

Note that it is only the long-term trend which counts. Short period fluctuations in weight over a few days are mostly the result of a loss of fluid.

The interpretation of the results will depend on changes in your:

overall weight and body fat, body water and muscle percentages, as well as on the period during which these changes take place. Rapid changes within days must be distinguished from medium term changes (over weeks) and long term changes (months).

A basic rule is that short term changes in weight almost exclusively represent changes in water content, whereas medium and long term changes may also involve the fat and muscle percentages.

- If your weight reduces over the short term, but your body fat percentage increases or remains the same, you have merely lost water – e.g. after a training session, sauna session or a diet restricted only to rapid weight loss.
- If your weight increases over the medium term and the body fat percentage falls or stays the same, then you could have built up valuable muscle mass.
- If your weight and body fat percentage fall simultaneously then your diet is working – you are losing fat mass. Ideally you should support your diet with physical activity, fitness or power training. By this means, you can increase your muscle percentage over the medium term.

Tips for using the scale

Important when measuring body fat/body water/muscle percentage:

- The measurement may only be carried out while barefoot and it is helpful if the soles of your feet are slightly damp. Completely dry soles can result in unsatisfactory results, as they have insufficient conductivity.
- Stand still during the measurement.
- Wait several hours (6-8) after unusually strenuous activity.
- Wait approx. 15 minutes after getting out of bed so that the water in your body can be distributed.

The measurement is not reliable for:

- Children under approx. 10 years of age;
- Professional sporters and bodybuilders;
- Pregnant women;
- Persons with fever, undergoing dialysis, with symptoms of edema or osteoporosis;
- Persons taking cardiovascular medication. Persons taking vasodilating or vasoconstricting medications;
- Persons with substantial anatomical deviations in the legs relative to their total height (leg length considerably shorter or longer than usual).

Incorrect measurement

If the scale detects an error during weighing, «Err» appears in the display.

Possible causes of errors:	Remedy:
The scale was not activated before stepping on. When you stand on the scale before "0.0" appeared in the display, it will not operate correctly.	Activate the scale the right way and repeat the measurement.
The maximum load-bearing capacity of 150 kg. (330 lbs) was exceeded.	Only weigh the maximum permissible weight.
The electrical resistance between the electrodes and the soles of your feet is too high (e.g. with heavily callused skin).	Repeat weighing barefoot. Slightly moisten the soles of your feet if necessary. Remove the calluses on the soles of your feet if necessary.
You were not standing still on the scale.	Stand still during the measurement.
Your body fat lies outside the measurable range (less than 5% or greater than 50%).	Repeat weighing barefoot and slightly moisten the soles of your feet if necessary.

Possible causes of errors:	Remedy:
Your body water lies outside the measurable range (less than 43% or greater than 75%).	Repeat weighing barefoot and slightly moisten the soles of your feet if necessary.
The percentage of muscle lies outside the measurable range (depending on age and sex).	Repeat weighing barefoot and slightly moisten the soles of your feet if necessary.

And remember: Physical exercise is the basis for a healthy body.

Convert from Kg (kilogram)/ St (Stone)/ Lb (Pounds)

This scale can display measurements in kilograms, Stone and pounds. You can set this with the switch on the bottom of the scale.

Kg (kilograms), St (Stone) or Lb (pounds).

C L E A N I N G & M A I N T E N A N C E

Clean the scale occasionally with a damp cloth. Do not use aggressive cleaning agents and never submerge the scale. **Never submerge the scale in water or other liquids.**

GUARANTEE & SERVICE

- This product is guaranteed for a period of 2 years from the date of purchase against mechanical and electrical defect.
- Guarantee can only be given if proof of purchase (e.g. till receipt) can be shown. Therefore keep it in a safe place. The guarantee period starts at date of purchase.
- Guarantee is not valid in case of:
 - normal wear;
 - damage through misuse or when the instructions were not followed correctly;
 - lack of maintenance;
 - not following the operation- and safety instructions;
 - incompetent montage or repairs made by yourself or others;
 - using accessories and spare parts which are not suitable for the appliance;for assessment by our technical department.
- The guarantee is not valid for:
 - batteries;for assessment by our technical department.

Return the appliance, including proof of purchase and description of the complaint, stamped to:

Inventum Holland BV
Generatorstraat 17
3903 LH VEENENDAAL
The Netherlands

- Pack it carefully (preferably in the original carton). Ensure unit is clean.
- Enclose your name and address and quote model number PW890BG on all correspondence.

When the guarantee is no longer valid or when the period of guarantee has passed, we will charge any costs that flow out of an assessment and/or repairing.

On request we will supply a statement before repairing. Without requesting a statement, we assume that you, by sending the appliance, automatically order to repair the appliance.

A guarantee period of two months is valid for repairs done by our technical department and the used spare parts, counting from the repair date. When a malfunction occurs within this period which directly flows out of the repairs and the used spare parts, the consumer has the right of repairs without any costs.

We reserve the right to change this product, packaging and documentation without notice.

Inventum Holland BV
Generatorstraat 17, 3903 LH , Veenendaal, Nederland
Tel.: 0318-542222, Fax.: 0318-542022, Mail: helpdesk@martex.nl

Inventum Belgium BVBA
Xavier de Cocklaan 66 bus 8, 9830, Deurle, België
Tel: 09 220 30 20, Fax.: 09 220 30 22, Mail: info@inventumbelgium.be

Website: www.inventum.eu