LANAFORM® MisterMasso●●

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INTRODUCTION

Lanaform® is committed to enabling you to benefit from the latest progress in comfort and body care. The relaxing, stimulating and reinvigorating effect of massagers on the muscles helps improve considerably your health and your figure. Just a few minutes a day is all it takes to relieve stress, tiredness and muscular pain. You should not, however, use your massager as a substitute for your doctor's advice. Persistent pain or other symptoms should be understood as indicators of a problem which may be more serious. Consult your doctor as soon as they appear before continuing to use your massager.

Description of the massager :

1. Adjustable speed control
For a soft to intense massage
2. OFF/M/MH switch:
OFF : to switch off the device
M : massage only
MH : massage with infrared
3. Ergonomic handle
4. Infrared massaging head
5. 2 X 2 attachments

PLEASE NOTE: do not activate infrared heat when using accessories (attachments)

IMPORTANT PRECAUTIONS to read carefully before use

As for any electrical appliance, a few fundamental safety precautions must be taken before using your "MISTER MASS".

- Please check that the voltage indicated on the massager plate matches the voltage in your home (220-230 volts)
- This device is designed for brief periods of use. It can operate without interruption for 15 minutes after which it is essential to switch it off and leave the motor to cool.
- Do not leave the massager unattended when it is operating or plugged in.
- Pregnant women, persons having had a thrombosis, persons experiencing body pain whose cause has not been diagnosed, and persons fitted with a heart pacemaker are warned not to use the "MISTER MASS" without medical advice.
- Particular care is needed when the massager is used by or in the vicinity of children or disabled people.

- Always keep the massager out of the reach of children.
- Never use this massager if the plug or electrical cable are damaged or are not working properly, or if it has been dropped, damaged, immersed in water or another liquid.
- As with any other electrical device, you are advised not to use the “MISTER MASS” during a thunderstorm.
- Do not use the massager close to a bath, sink or other areas where it could come into contact with water or damp.
- Never immerse the massager in water, in your bath or in the shower.
- Keep the electrical cable well clear of any heating or cutting device.
- Unplug before cleaning.
- Check the massager’s temperature before using on a sensitive area.
- Do not use the “MISTER MASS” around the eyes or close to any sensitive part of the body or close to wounds or infected areas.
- Do not use the massager while you are sleeping or without continuously attending to it.
- Do not use the massager if aerosols are being used or oxygen administered in the vicinity.
- Do not use the massager under a cushion or pillow as this could cause it to overheat, resulting in a fire, electric shock or injury to persons.
- Never carry the massager by its electrical cable and do not use the cable as a handle.
- Before unplugging the massager, ensure that all the controls are in the OFF position.
- Only use the “MISTER MASS” for the purposes for which it is designed and as stated in this manual.

Instructions for use

1. The massager must be switched OFF (switch down) when you plug it in.
2. To start the “MISTER MASS”, use the “OFF/M/MH” switch located on the top of the handle. To adjust the massage intensity, turn the speed control to the desired position.
3. Move the massage head slowly across the surface to be massaged. Massage the surface for a brief period and then move the massager to the next area to be massaged. Never massage the same spot for more than 3 minutes.
4. You are advised not to massage for more than 15 minutes as that can cause the device to overheat.
5. When you have finished massaging, switch “OFF”, unplug and store in a dry place.

TWO-YEAR UNLIMITED WARRANTY

LANAFORM® guarantees the “MISTER MASS” against all material or manufacturing defects for a period of two years from the date of purchase, provided it is used for normal domestic purposes and subject to the exceptions as stated hereunder.

This LANAFORM® warranty does not cover damage caused by improper use, normal wear and tear, accident, the fitting of any unauthorised accessory, deterioration of the product, careless handling during transportation or any other cause beyond the control of LANAFORM®

All warrantees, including but not limited to warrantees of correctness and fair and merchantable quality, are limited for a maximum of two years from the date of purchase of the product.

On reception of a defective product under warranty, LANAFORM® will repair or if applicable replace the massager and return it to you. The warranty is exercised by LANAFORM® Service Center only. The repair of this product by any other person than Lanaform® Service Center will result in the cancellation of the warranty.

Address :

SA LANAFORM NV
Zoning de Cornémont
B-4141 LOUVEIGNE – BELGIUM
Tel 32/4 360 92 91 -- FAX 32/4 360 97 23
info@lanaform.com

MAINTENANCE:

1. **STORAGE:** Place the massager in its box in a safe, dry and cool place. Do not allow the device to come into contact with cutting or pointed objects which could damage the external surface. Do not wind the electrical cable around the massager. (This could cause the cable to become worn and even break).
2. **CLEANING:** Unplug the massager and leave to cool before cleaning.
3. **DO NOT IMMERSE THE "MISTER MASS" IN A LIQUID TO CLEAN IT.** Use a mild detergent. Dry by dabbing with a clean, dry cloth.
4. **NEVER USE** abrasive detergents, a brush, petrol, kerosene, cleaning products for glass or furniture, paint thinner or water to clean the massager.

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