SALTER

PERSONAL SCALES

Instructions and Guarantee

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PERSONENWAAGE

BALANZA PERSONAL

BILANCIA PESAPERSONA

BILANÇÀ PESAPERSONA

BADEROMSVEKT

PERSONENWEEGSCHAAL

HENKILÖVAAAT

PERSONVÅG

PERSONVÆGT

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OSOBNÍ VÁHA

KİŞİSEL TERAZİ

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OSOBNÉ VÁHY

ULTIMATE ACCURACY PROGRESS TRACKER

This Salter Ultimate Accuracy Progress Tracker features 50g increments for precise weight management. The Progress Tracker aids effective weight management. It displays your current weight, start weight and the difference, indicating total weight change. It also displays your previous weight reading and the difference between current weight indicating recent progress.

PREPARING YOUR SCALE

1. Open the battery compartment on the scale underside.
2. Insert batteries (1 x CR2032) observing the polarity signs (+ and -) inside the battery compartment.
3. Close the battery compartment.
4. Select kg, st or lb weight mode by the switch on the underside of the scale.
5. For use on carpet attach enclosed carpet feet.

WEIGHT READING ONLY

1. Position scale on a firm flat surface.
2. Tap the platform centre (vibration activates your scale) and remove your foot.
3. Wait until zero is displayed.
4. Step on and stand very still while the scale computes your weight.
5. Your weight is displayed for a few seconds then the scale will switch off.

SET START WEIGHT

1. Position scale on a firm flat surface.
2. Select user number 1-4
3. Press and hold the user number for 3 seconds. "Start Weight 0.00" will be displayed then the scale will turn off.
4. Select user number 1-4.
5. Wait until zero is displayed.
6. Step on and stand very still while the scale computes your weight.
7. The current "Weight" reading will be set as "Start Weight" automatically the next time you use the scale.
8. Repeat the procedure for a second user or to update the "Start Weight" reading.

PROGRESS TRACKING

1. Position scale on a firm flat surface.
2. Select user number 1-4
3. Wait until zero is displayed.
4. Step on and stand very still while the scale computes your weight.
5. Your current "Weight" is displayed for a few seconds followed by "Start Weight", "Start Weight Difference", "Previous Weight" and "Previous Weight Difference".

Note: A down arrow indicates weight less than "Start Weight"/"Previous Weight". An up arrow indicates weight more than "Start Weight"/"Previous Weight".

6. Your current "Weight" is displayed a second time then the scale will turn off.

WARNING INDICATORS

- Err** Unstable weight. Stand still.
O-Ld Weight exceeds maximum capacity.
Lo Replace battery.

ADVICE FOR USE AND CARE

- Always weigh yourself on the same scale placed on the same floor surface. Do not compare weight readings from one scale to another as some differences will exist due to manufacturing tolerances.
- Placing your scale on a hard, even floor will ensure the greatest accuracy and repeatability.
- Weigh yourself at the same time each day, before meals and without footwear. First thing in the morning is a good time.
- Your scale rounds up or down to the nearest increment. If you weigh yourself twice and get two different readings, your weight lies between the two.
- Clean your scale with a damp cloth. Do not use chemical cleaning agents.
- Do not allow your scale to become saturated with water as this can damage the electronics.
- Treat your scale with care - it is a precision instrument. Do not drop it or jump on it.

TROUBLESHOOTING GUIDE

If you experience any difficulty in the use of your scale:

- Check the battery is correctly fitted.
- Check that you have selected your preference of stones/pounds, kilograms or pounds.
- Check that the scale is on a flat, level floor and not touching against a wall.
- If when you use the scale there is no display at all or if 'Lo' is displayed - try a new battery.
- If 'O-Ld' is displayed the scale has been overloaded.

WEEE EXPLANATION

 This marking indicates that this product should not be disposed with other household wastes throughout the EU. To prevent possible harm to the environment or human health from uncontrolled waste disposal, recycle it responsibly to promote the sustainable reuse of material resources. To return your used device, please use the return and collection systems or contact the retailer where the product was purchased. They can take this product for environmental safe recycling.

BATTERY DIRECTIVE

 This symbol indicates that batteries must not be disposed of in the domestic waste as they contain substances which can be damaging to the environment and health. Please dispose of batteries in designated collection points.

GUARANTEE

This product is intended for domestic use only. Salter will repair or replace the product, or any part of this product, free of charge if within 15 years of the date of purchase, it can be shown to have failed through defective workmanship or materials. This guarantee covers working parts that affect the function of the scale. It does not cover cosmetic deterioration caused by fair wear and tear or damage caused by accident or misuse. Opening or taking apart the scale or its components will void the guarantee. Claims under guarantee must be supported by proof of purchase and be returned carriage paid to Salter (or local Salter appointed agent if outside the UK). Care should be taken in packing the scale so that it is not damaged while in transit. This undertaking is in addition to a consumer's statutory rights and does not affect those rights in any way. For UK Sales and Service contact HoMedics Group Ltd, PO Box 460, Tonbridge, Kent, TN9 9EW, UK. Helpline Tel No: (01732) 360783. For Ireland, please contact Petra Brand Masters, Unit J4 Maynooth Business Campus, Maynooth, Co. Kildare, Ireland. Tel +00 353 (0) 1 6510660. e-mail sales@petrabrandmasters.ie.
www.salterhousewares.com/servicecentres

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1. Position scale on a firm flat surface.
2. Select user number 1-4
3. Wait until zero is displayed.
4. Step on and stand very still while the scale computes your weight.
5. Your current "Weight" is displayed for a few seconds followed by "Start Weight", "Start Weight Difference", "Previous Weight" and "Previous Weight Difference".
Note: A down arrow indicates weight less than "Start Weight"/"Previous Weight". An up arrow indicates weight more than "Start Weight"/"Previous Weight".
6. Your current "Weight" is displayed a second time then the scale will turn off.

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