SANITAS

SBF 11

(D) Gebrauchsanleitung

Diagnosewaage

(GB) Operating Instructions

Diagnostic scale

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Bilancia diagnostica

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Please read these instructions for use carefully and keep them for later use, be sure to make them accessible to other users and observe the information they contain.

1. Signs and symbols

The following symbols appear in these instructions.

Warning

Warning instruction indicating a risk of injury or damage to health.

Important

Safety note indicating possible damage to the unit/accessory.

Note

Note on important information.

2. Safety notes

Please read these instructions for use carefully and keep them for later use, be sure to make them accessible to other users and observe the information they contain.

WARNING

- **The scales must not be used by persons with medical implants (e.g. heart pacemakers) Otherwise their function could be impaired.**
- Do not use during pregnancy.
- Do not stand on the outer edge to one side of the scale otherwise it may tip!
- Keep children away from packaging materials (risk of suffocation).
- Caution! Do not step onto the scale with wet feet or when the scale's surface is damp - danger of slipping!

Notes on handling batteries

- Swallowing batteries can be extremely dangerous. Keep the batteries and scale out of the reach of small children. Should a child swallow a battery, seek medical assistance immediately.
- Replace weak batteries before they discharge completely.
- Batteries must not be recharged, taken apart, thrown into an open fire or short circuited.
- Leaking batteries may damage the device. If you do not intend to use the device for longer periods, remove the batteries from the battery compartment.
- If a battery has leaked, put on protective gloves and clean the battery compartment with a dry cloth.
- Batteries can contain toxins that are harmful to health and the environment. Always dispose of batteries in accordance with applicable legal regulations. Do not dispose of batteries with the normal household waste.

General notes

- The unit is for personal use only and is not intended for medical or commercial applications.
- Note that technical tolerances mean that results may vary, because the scale is not calibrated for professional medical use.
- The scale can accept a maximum load of 150 kg (330 lb, 23 st/8 lb). The weight measurement is displayed in increments of 100 g (0.2 lb).
- The measuring results of the body fat, body water and muscle percentage are shown in 0.1% increments.
- On delivery, the scale is set to „cm“ and „kg“. At the back of the scales is a switch to change between „pounds“ and „stones“.
- Place the scale on a firm level floor; a firm floor covering is required for correct measurement.
- Repairs may only be carried out by customer service or authorised dealers. Before submitting any complaint, first check the batteries and replace them if necessary.

- Should you have any questions concerning the use of our devices, please contact your dealer or customer service.
- All scales comply with EC Directive 2004/108/EC. Should you have any questions concerning the use of our devices, please contact your dealer or customer service.

Storage and maintenance

The precision of the measured values and service life of the unit depend on careful use:

IMPORTANT

- The scale should be cleaned occasionally with a damp cloth. Do not use abrasive detergents and never immerse the scale in water.
- Ensure that no liquid gets on the scale. Never dip the scale into water. Never wash it under running water.
- Do not place any objects on the scale when it is not being used.
- Protect the scale from impacts, moisture, dust, chemicals, large changes in temperature and sources of heat that are too close (oven, radiator).
- Do not press the buttons using excessive force or with pointed objects.
- Do not expose the scale to high temperatures or strong electromagnetic fields (e.g. mobile phones).

Batteries

WARNING

- Do not throw batteries into fire. Explosion hazard!

Repairs


IMPORTANT

- Never open or attempt to repair the unit yourself, as otherwise proper function is no longer guaranteed. Failure to observe this regulation shall void the warranty.
- If you need to have the unit repaired, contact customer service or an authorised dealer.

3. Unit description

Overview

1. Electrodes
2. Display
3. Up key
4. ON button
5. SET button
6. Down key

4. Initial use

Inserting batteries

If provided, pull out the battery insulating strip at the battery compartment or remove the protective film from the battery and insert the battery with the correct polarity. If the scale fails to operate, remove the battery completely and reinsert it.

Your scale is provided with a "Change battery" indicator. If the scale is operating with a battery which is too weak, "LO" appears on the display and the scale will switch off automatically. In this case, the battery must be exchanged (1 x3 V Lithium CR 2032).

Setting up the scale

Set up the scale on a level and solid surface. A solid surface is a prerequisite for a correct measurement.

5. Information

The measuring principle

This scale operates according to the B.I.A. principle (bioelectric impedance analysis). This enables the measurement of physical relationships within seconds by means of an undetectable, completely harmless electric current. The body fat percentage and other physical relationships in the body can be determined by measuring the electrical resistance (impedance) and calculating constants and individual parameters such as age, height, gender and degree of activity.

Muscle tissue and water have good electrical conductivity, and therefore lower resistance. Bones and fatty tissue, on the other hand, have low conductivity as fat cells and bones conduct hardly any current as a result of their very high resistance.

Please be aware that values obtained from the diagnostic scale represent only an approximation of actual analytical medical data. Only a specialist physician can accurately determine body fat, body water, muscle percentage and bone structure using medical procedures (e.g. computed tomography).

General tips

- If possible, always weigh yourself at the same time of day (preferably in the morning), after using the toilet, on an empty stomach and unclothed, in order to achieve comparable results.
- Important for the measurement: only measure body fat when barefoot and with the soles of the feet slightly moist.
Completely dry soles can result in unsatisfactory measurements due to inadequate conductivity.
- Stand still and upright during the measurement.
- Wait a few hours after unaccustomed physical activity.
- Wait approx. 15 minutes after getting up in the morning to allow the water in the body to distribute.
- Remember that only long-term trends are important. Short-term changes in weight within a few days are usually caused by a loss of fluid. Body water plays an important role towards our general wellbeing.

Limitations

When measuring body fat and other values, deviating and implausible results may occur in:

- children under approx. 10 years,
- competitive athletes and body builders,
- pregnant women,
- persons with fever, undergoing dialysis treatment or with symptoms of edema or osteoporosis,
- persons taking cardiovascular medication (affecting the heart and vascular system),
- persons taking vascular dilating or vascular constricting medication,
- persons with considerable anatomic deviations of the legs in relation to total body size (length of the legs considerably shortened or lengthened).

6. Settings

Enter your personal information before using the scale.

On delivery, the scale is set to „cm“ and „kg“. On the rear of the scale is a switch that you can use to change to 'pounds/stones'.

Setting user data

You must enter personal user parameters in order to measure your body fat percentage and other physical data.

The scale has 10 user memory locations which allows, for example, you and your family members to save personal settings and retrieve them later.

- Press 'ON'.
- Wait until '0.0' is displayed.
- Press 'SET'.
- Select the required memory space with the ▼/▲ key and confirm with the 'SET' key
- Now you can adjust the following settings:

User data	Parameters
Age	10 to 100 Jahre
Body height	100 to 220 cm (3' 03" to 7' 03")
Sex	man (♂), woman (♀)
Degree of activity	1 to 5

- Modifying values: Press ▲ or ▼ and hold for more rapid progress.
- Confirming entries: Press the 'SET' button.
- After the values are set, '0.0' appears in the display as well as the entered user data.
- Now the scale is ready to measure your weight. Now step onto the scale. Stand still on the scale with your weight distributed evenly between both legs. The scale immediately begins to measure your weight.

If no further action is taken, the scale switches itself off.

Degrees of activity

Selection of the degree of activity must refer to the medium and long term.

Degree of activity	Physical activity
1	None.
2	Low: A small amount of light physical effort (e.g. short walks, light garden work, gymnastic exercises).
3	Medium: Physical effort for 30 minutes at least 2 to 4 times a week.
4	High: Physical effort for 30 minutes at least 4 to 6 times a week.
5	Very high: Intensive physical effort, intensive training or hard physical work for at least one hour daily.

All parameters must be entered before the body fat and other values can be entered.

7. Operation

Place the scale on a firm level floor (no carpet); a firm floor covering is required for correct measurement.

Only weight measurement

Press the 'ON' button. In the display, '8888' appears (fig. 1), then '0.0' (fig. 2).

Now the scale is ready to measure your weight. Now step onto the scale. Stand still on the scale with your weight distributed evenly between both legs. The scale immediately begins to measure your weight.

fig. 1

Weighing, carrying out diagnostics

- Press the 'ON' button. In the display, '8888' appears (fig. 1), then '0.0' (fig. 2).
- Press the ▲ or ▼ buttons several times to select the memory location in which your basic personal data is saved and confirm by pressing 'SET'. The last measuring results are displayed in sequence until '0.0' appears.
- If no measurements have been saved, '----' is displayed.
- Step onto the scale barefoot and make sure that you are standing still with your weight evenly distributed between both legs on the stainless steel electrodes.

fig. 2

fig. 3

Note: Your feet, legs, calves and thighs must not touch each other. Otherwise the measurement cannot be performed correctly. The scale immediately begins to measure your weight.

fig. 4

The scale begins to take measurements straight away. First, the weight is displayed. Whilst the other parameters are being measured, an 'arrow' at the bottom edge of the LCD' appears (fig. 3).

Your measured weight is displayed shortly after this (Fig 4).

The following appears:

1. Display:

- Weight in kg (lb, st)
- Body fat in % with interpretation BF
- Body water in % $\approx$

2. Display:

- Muscle percentage in % $\text{—} \text{■} \text{—}$
- Bone mass in kg $\text{—} \text{■} \text{—}$
- Active metabolic rate in kcal

Switching off the scale

The scale automatically switches off after approx. 10 seconds.

Press the 'ON' button to switch off manually.

8. Evaluation of results

Body fat percentage

Body fat content is evaluated at the bottom of the display. The following body fat percentages are guide values (contact your physician for further information).

Man

				
Age	low	normal	high	very high
10-14	<11 %	11-16 %	16,1-21 %	>21,1 %
15-19	<12 %	12-17 %	17,1-22 %	>22,1 %
20-29	<13 %	13-18 %	18,1-23 %	>23,1 %
30-39	<14 %	14-19 %	19,1-24 %	>24,1 %
40-49	<15 %	15-20 %	20,1-25 %	>25,1 %
50-59	<16 %	16-21 %	21,1-26 %	>26,1 %
60-69	<17 %	17-22 %	22,1-27 %	>27,1 %
70-100	<18 %	18-23 %	23,1-28 %	>28,1 %

Woman

				
Age	low	normal	high	very high
10-14	<16 %	16-21 %	21,1-26 %	>26,1 %
15-19	<17 %	17-22 %	22,1-27 %	>27,1 %
20-29	<18 %	18-23 %	23,1-28 %	>28,1 %
30-39	<19 %	19-24 %	24,1-29 %	>29,1 %
40-49	<20 %	20-25 %	25,1-30 %	>30,1 %
50-59	<21 %	21-26 %	26,1-31 %	>31,1 %
60-69	<22 %	22-27 %	27,1-32 %	>32,1 %
70-100	<23 %	23-28 %	28,1-33 %	>33,1 %

A lower value is often found in athletes. Depending on the type of sports, training intensity and physical constitution, values may result which are below the recommended values stated. It should, however, be noted that there could be a danger to health in the case of extremely low values.

Body water percentage

The body water percentage is normally within the following ranges:

Man

Age	poor	gooe	very good
10-100	<50 %	50-65 %	>65 %

Woman

Age	poor	gooe	very good
10-100	<45 %	45-60 %	>60 %

Body fat contains relatively little water. Therefore persons with a high body fat percentage have body water percentages below the recommended values. With endurance athletes, however, the recommended values could be exceeded due to low fat percentages and high muscle percentage.

This scale is unsuitable for measuring body water in order to draw medical conclusions concerning age-related water retention, for example. If necessary ask your physician. Basically, a high body water percentage should be the aim.

Muscle percentage

The muscle percentage is normally within the following ranges:

Man

Age	low	normal	high
10-14	<44 %	44-57 %	>57 %
15-19	<43 %	43-56 %	>56 %
20-29	<42 %	42-54 %	>54 %
30-39	<41 %	41-52 %	>52 %
40-49	<40 %	40-50 %	>50 %
50-59	<39 %	39-48 %	>48 %
60-69	<38 %	38-47 %	>47 %
70-100	<37 %	37-46 %	>46 %

Woman

Age	low	normal	high
10-14	<36 %	36-43 %	>43 %
15-19	<35 %	35-41 %	>41 %
20-29	<34 %	34-39 %	>39 %
30-39	<33 %	33-38 %	>38 %
40-49	<31 %	31-36 %	>36 %
50-59	<29 %	29-34 %	>34 %
60-69	<28 %	28-33 %	>33 %
70-100	<27 %	27-32 %	>32 %

Bone mass

Like the rest of our body, our bones are subject to the natural development, degeneration and ageing processes.

Bone mass increases rapidly in childhood and reaches its maximum between 30 and 40 years of age. Bone mass reduces slightly with increasing age. You can reduce this degeneration somewhat with healthy nutrition (particularly calcium and vitamin D) and regular exercise. With appropriate muscle building, you can also strengthen your bone structure. Note that this scale will not show you the calcium content of your bones, but will measure the weight of all bone constituents (organic substances, inorganic substances and water). Little influence can be exerted on bone mass, but it will vary slightly within the influencing factors (weight, height, age, gender). no recognised guidelines or recommendations relating to bone mass measurement.

ATTENTION:

Please do not confuse bone mass with bone density.

Bone density can be determined only by means of a medical examination (e.g. computer tomography, ultrasound). It is therefore not possible to draw conclusions concerning changes to the bones and bone hardness (e.g. osteoporosis) using this scale.

AMR

The active metabolic rate (AMR) is the amount of energy required daily by the body in its active state. The energy consumption of a human being rises with increasing physical activity and is measured on the diagnostic scale in relation to the degree of activity entered (1-5).

To maintain your existing weight, the amount of energy used must be reintroduced into the body in the form of food and drink. If less energy is introduced than is used over a longer period of time, your body will obtain the difference largely from the amount of fat stored and your weight will decrease. If, on the other hand, over a longer period of time more energy is introduced than the total active metabolic rate (AMR) calculated, your body will be unable to burn off the excess energy, and the excess will be stored in the body as fat and your weight will increase.

Results in relation to time

-  Remember that only long-term trends are important. Short-term fluctuations in weight over a few days are usually the result of a loss of fluid.

The interpretation of the results will depend on changes in your overall weight and body fat, body water and muscle percentages, as well as on the period during which these changes take place. Sudden changes within days must be distinguished from medium term changes (over weeks) and long term changes (months).

A basic rule is that short term changes in weight almost exclusively represent changes in water content, whereas medium and long term changes may also involve the fat and muscle percentages.

- If your weight reduces over the short term, but your body fat percentage increases or remains the same, you have merely lost water – e.g. after a training session, sauna session or a diet restricted only to rapid weight loss.
- If your weight increases over the medium term and the body fat percentage falls or stays the same, then you could have built up valuable muscle mass.
- If your weight and body fat percentage fall simultaneously then your diet is working – you are losing fat mass.
- Ideally you should support your diet with physical activity, fitness or power training. This way you can increase your muscle percentage over the medium term.
- Body fat, body water or muscle percentages should not be added (certain elements of muscle tissue also contain body water).

9. Disposal

Standard and rechargeable batteries should not be disposed of separately from the household waste. As a consumer, you are legally obliged to return used batteries for proper disposal. You can hand in your used batteries at public collection points in your district or sales outlets where batteries of this type are sold.


NOTE:

The codes below are printed on batteries containing harmful substances:

Pb = Battery contains lead,

Cd = Battery contains cadmium,

Hg = Battery contains mercury.

For environmental reasons, do not dispose of the scale with the batteries in the household waste at the end of its useful life. Dispose of the unit at a suitable local collection or recycling point. Observe the local regulations for material disposal.

Dispose of the unit in accordance with EC Directive 2002/96/EC – WEEE (Waste Electrical and Electronic Equipment).

If you have any questions, please contact the local authorities responsible for waste disposal.

10. What if there are problems?

If the scale detects an error when measuring, the following is displayed.

Display	Cause	Remedy
0-Ld	The maximum load-bearing capacity of 150 kg was exceeded.	The load must not exceed 150 kg.
LO	The scale batteries are almost flat.	Replace the batteries of scale.
Err	The fat percentage is outside the measurable range (less than 3% or greater than 55%). The water percentage is outside the measurable range (less than 25% or greater than 75%).	Repeat the measurement barefoot or, if necessary, moisten the soles of your feet slightly.

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