SANITAS

SGS 43

- (D) **Gebrauchsanweisung**
Glaswaage
- (GB) **Instruction for use**
Glass scale
- (F) **Mode d'emploi**
Pèse-personne en verre
- (E) **Instrucciones para el uso**
Báscula de vidrio
- (RUS) **Инструкция по применению**
Стеклянные весы

Service-Adresse:
Hans Dinslage GmbH
Riedlinger Str. 28
88524 Uttenweiler
GERMANY
Tel.: +49 (0) 73 74 - 91 57 66
Fax: +49 (0) 73 74 - 92 07 23
E-mail: service@sanitas-online.de

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Contents

1. Getting to know your instrument	10
2. Signs and symbols.....	10
3. Safety notes.....	11
4. Unit description.....	12
5. Initial use	12
6. Setting up the control panel.....	13
7. Information.....	13
8. Operation.....	13
9. Replacing batteries	14
10. Cleaning and care of the unit.....	14
11. Disposal	15
12. What if there are problems?	15
13. Technical specifications	16

Included in delivery

- Glass scale
- Removable control unit with display
- Wall mount
- 4 x 1.5 V AAA batteries (scale)
- 2 x CR2032 lithium batteries (control unit)
- These instructions for use

Please read these instructions for use carefully and keep them for later use, be sure to make them accessible to other users and observe the information they contain. This diagnostic scale is designed to make a significant contribution to your health.

1. Getting to know your instrument

Function of the unit

The digital scale is used for measuring weights. The scale has a display that receives the measured values by infrared transmission. Weight is displayed in 100 g graduations. The scale is designed for private use.

This scale also has the following functions:

- Switch between kilograms, pounds („LB“) and stone („ST“).
- Automatic shutoff function
- Battery change indicator for weak batteries

2. Signs and symbols

The following symbols appear in these instructions.

Warning

Warning instruction indicating a risk of injury or damage to health.

Important

Safety note indicating possible damage to the unit/accessory.

Note

Note on important information.

3. Safety notes

Please read these instructions for use carefully and keep them for later use, be sure to make them accessible to other users and observe the information they contain.

- Do not stand on the outer edge to one side of the scale otherwise it may tip!
- Keep children away from packaging materials (risk of suffocation).

Notes on handling batteries

- Swallowing batteries can be extremely dangerous. Keep the batteries and scale out of the reach of small children. Should a child swallow a battery, seek medical assistance immediately.
- Replace weak batteries before they discharge completely.
- Always replace all batteries at the same time and use batteries of the same type.
- Batteries must not be recharged, taken apart, thrown into an open fire or short circuited.
- Leaking batteries may damage the device. If you do not intend to use the device for longer periods, remove the batteries from the battery compartment.
- If a battery has leaked, put on protective gloves and clean the battery compartment with a dry cloth.
- Batteries can contain toxins that are harmful to health and the environment. Always dispose of batteries in accordance with applicable legal regulations. Do not dispose of batteries with the normal household waste.

General notes

- The unit is for personal use only and is not intended for medical or commercial applications.
- Caution! Do not step onto the scale with wet feet or when the scale's surface is damp - danger of slipping!
- Note that technical tolerances mean that results may vary, because the scale is not calibrated for professional medical use.
- The scale can accept a maximum load of 150 kg (330 lb, 23 st). The weight measurement is displayed in increments of 100 g (0.2 lb).
- Place the scale on a firm level floor; a firm floor covering is required for correct measurement.
- Protect the unit against hard knocks, moisture, dust, chemicals, major temperature fluctuations and heat sources which are too close (stoves, heating radiators).
- Repairs may only be carried out by customer service or authorised dealers. Before submitting any complaint, first check the batteries and replace them if necessary.
- Should you have any questions concerning the use of our devices, please contact your dealer or customer service.
- All scales comply with EC Directive 2004/108/EC. Should you have any questions concerning the use of our devices, please contact your dealer or customer service.

Storage and maintenance

The precision of the measured values and service life of the unit depend on careful use:

Important

- The unit should be cleaned occasionally. Do not use abrasive detergents and never immerse the unit in water.
- Ensure that no liquid gets on the scale. Never dip the scale into water. Never wash it under running water.
- Do not place any objects on the scale when it is not being used.
- Protect the scale from impacts, moisture, dust, chemicals, large changes in temperature and sources of heat that are too close (oven, radiator). Do not press the buttons using excessive force or with pointed objects.
- Do not expose the scale to high temperatures or strong electromagnetic fields (e.g. mobile phones).

Batteries

Warning

- Swallowing batteries can be extremely dangerous. Keep the batteries out of the reach of small children. If a battery is swallowed, get immediate medical help.
- Do not throw batteries into fire. Explosion hazard!

Important

- Batteries must not be recharged or reactivated by other means, taken apart or short-circuited.
- Remove batteries if the unit will not be used for an extended period of time. This prevents damage that may occur from leaks.

Repairs


Important

- Never open or attempt to repair the unit yourself, as otherwise proper function is no longer guaranteed. Failure to observe this regulation shall void the warranty.
- If you need to have the unit repaired, contact customer service or an authorised dealer.

4. Unit description

Overview

- 1 Control panel with display
- 2 ON/OFF button
- 3 Receiver
- 4 Scale
- 5 Transmitter
- 6 LED
- 7 Wall mount


5. Initial use

Inserting batteries

If present, pull the battery insulating strip off the battery compartment cover or remove the battery's protective film and insert the battery in the scale and the control panel with the correct polarity. Note the graphic in the battery compartment.

Control unit


2 x CR2032 lithium batteries

Scale

4 x 1.5 V AAA batteries

If the scale fails to operate, remove the batteries completely and reinsert them.

For instructions on changing the batteries, see p. 14.

Changing the weight unit

On delivery the scale is set to the unit „kg“. A switch located on the back of the control panel enables you to change to „pounds“ (lb) and „stone“ (st).

Setting up the scale

Set up the scale on a level and solid surface. A solid surface is a prerequisite for a correct measurement.

6. Setting up the control panel


Positioning the control panel

The control panel can be held in the hand, attached to the wall mount or be left on the scale.

Wall mounting

The wall mount can be attached to the wall with the double-sided adhesive strip. The surface of the wall must be clean, plane and grease-free. Remove the foil from the adhesive strip on the back of the wall mount. Press the wall mount firmly against the wall.

- Make sure that the wall mount is straight.
- Note that the control panel and the scale must also be aligned to each other.

7. Information

General tips




- If possible, always weigh yourself at the same time of day (preferably in the morning), after using the toilet, on an empty stomach and unclothed, in order to achieve comparable results.
- Stand still and upright during the measurement.

8. Operation

Weighing

Place the scale on a firm level floor (no carpet); a firm floor covering is required for correct measurement.

Make sure that the control panel and the scale are in line of sight to each other. Do not cover the receiver and transmitter. The maximum distance between control panel and scale is 2 m.

Option 1: switch on control panel, then scale platform

Press the „ON/OFF“ button on the display module. The display shows „8888“ (Fig. 1).

Audibly, briefly and firmly tap the scale with your foot to activate the vibration sensor. The LED first shows red and then switches to green after about 5 seconds. The display simultaneously shows „0.0“ (Fig. 2).

Now the scale is ready to measure your weight. Now step onto the scale. Stand still on the scale with your weight distributed evenly between both legs. The scale immediately begins to measure your weight. The display shows „---“ (Fig. 3).

Fig. 1

Fig. 2

Your measured weight is displayed shortly after this (Fig 4).

The scale switches off 10 seconds after you step down from it. The result remains visible during this time.

Make sure that you always switch on the scale first and wait for the display to show „0.0“ (Fig. 2) before stepping on it.

Fig. 3

Fig. 4

Option 2: switch on scale platform, then control panel

Audibly, briefly and firmly tap the scale with your foot to activate the vibration sensor. The LED first shows red and then switches to green after about 5 seconds.

Press the „ON/OFF“ button on the display module. If the control panel is receiving the signal from the scale it shows „0.0“ (Fig. 2).

Now the scale is ready to measure your weight. Now step onto the scale. Stand still on the scale with your weight distributed evenly between both legs. The scale immediately begins to measure your weight. The display shows „----“ (Fig. 3).

Your measured weight is displayed shortly after this (Fig 4).

The scale switches off 10 seconds after you step down from it. The result remains visible during this time.

Make sure that you always switch on the scale first and wait for the display to show „0.0“ (Fig. 2) before stepping on it.

Switching off the scale

The scale and the control panel always switch off automatically after 10 seconds.

Press the „ON/OFF“ button to switch off the display module manually.

9. Replacing batteries

The control panel is equipped with a „replace battery“ indicator.

The display shows „L0-1“ and „L0-2“ when the batteries are beginning to run down and the scale switches off automatically. The batteries must be replaced as soon as possible. If the display shows „L0-1“, the batteries in the scale must be replaced (4 x 1.5 V AAA). If the display shows „L0-2“, the batteries in the control panel must be replaced (2 x CR2032 lithium batteries).

The batteries in the scale must also be replaced if the LED on the scale does not light after activating the vibration sensor (4 x 1.5 V AAA batteries).

Note:

- When changing batteries, always use batteries of the same type, brand and capacity.
- Do not use rechargeable batteries.
- Use batteries that are free of heavy metals.

10. Cleaning and care of the unit

The unit should be cleaned occasionally.

Clean using a damp cloth, to which you can apply a little detergent if necessary. Do not use abrasive detergents and never immerse the unit in water.

Important

- Do not use aggressive solvents or cleaning agents!
- Never immerse the unit in water!
- Do not wash the unit in a dishwasher!

11. Disposal

Standard and rechargeable batteries should not be disposed of separately from the household waste. As a consumer, you are legally obliged to return used batteries for proper disposal. You can hand in your used batteries at public collection points in your district or sales outlets where batteries of this type are sold.


Note:

The codes below are printed on batteries containing harmful substances:

Pb = Battery contains lead,

Cd = Battery contains cadmium,

Hg = Battery contains mercury.

For environmental reasons, do not dispose of the scale with the batteries in the household waste at the end of its useful life. Dispose of the unit at a suitable local collection or recycling point.

Observe the local regulations for material disposal.

Dispose of the unit in accordance with EC Directive **2002/96/EC – WEEE** (Waste Electrical and Electronic Equipment).

If you have any questions, please contact the local authorities responsible for waste disposal.

12. What if there are problems?

If the scale detects an error when measuring, the following is displayed.

Display	Cause	Remedy
0-Ld	The maximum load-bearing capacity of 150 kg was exceeded.	The load must not exceed 150 kg.
Lo-1	The scale batteries are almost flat.	Replace the batteries of scale.
Lo-2	The control panel batteries are almost flat.	Replace the batteries of control panel.
Err	Faulty or no transmission between scale platform and control panel.	Place the scale on a flat surface and activate it correctly (wait for the display to show „0.0“) and repeat the measurement (see page 14).
No display	The batteries in the control panel are entirely flat.	Replace the batteries (see page 13).
	The batteries in the control panel are not correctly installed.	Check that the polarity is correct (see page 12).

No measurement possible

Possible errors	Remedy
The scale was not activated before being stepped on. If you step onto the scale before „0.0“ appears in the display, it will not function correctly.	Activate the scale correctly (wait for „0.0“) and repeat the measurement. Switch on the control panel and the scale before the measurement. If necessary, wait 30 seconds and then repeat the measurement.
Interference between scale and control panel.	Remove and replace the batteries in the scale and control panel or make sure that the transmitter and receiver are not covered.
The control panel is too far from the scale.	Reduce the distance between the scale and control panel.
The scale batteries are flat.	Replace the batteries (see page 14).
The connection between control panel and scale is interrupted.	Do not cover the receiver and transmitter.

Incorrect weight measurement

Possible errors	Remedy
Scale is on carpet.	Place scale on a flat, solid surface.
Scale zero setting is incorrect.	Switch on scale and wait until it switches off automatically (LED off). Then repeat measurement.

13. Technical specifications

Scale batteries	4 x 1.5 V AAA
Control unit batteries	2 x CR2032 lithium batteries
Measuring range	1 kg to 150 kg
Weight	100 g
Transmission	Infrared transmission

Subject to technical modifications.

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