# SOEHNLE

Body  
Balance  
Adria

Bedienungsanleitung  
Operating Instructions  
Notice d'utilisation  
Istruzioni per l'uso  
Gebruiksaanwijzing  
Instrucciones manejo

Körperwasser- & Körperfett-Analyse-Waage • Body water and body fat monitoring scale • Balance d'analyse de l'eau et des graisses du corps • Bilancia per l'analisi del grasso e dell'acqua nel corpo • Lichaamswater & lichaamsvet analyseweegschaal • Báscula de análisis de líquido y grasa corporal

**Inhalt****Body Balance Adria von SOEHNLE**

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## Body Balance Adria by SOEHNLE

By choosing the Body Balance Adria, you are making an important contribution to your health care. This electronic body monitoring scale performs the following functions:

- ▶ Precise body water, body fat and weight measurement
- ▶ individual interpretation of results with personal advisor function (see page 17)
- ▶ 4-person memory for individual specific data storage

### Important advice

We recommend that you consult a doctor if you are overweight and are trying to reduce your weight, or if you are underweight and are trying to put on weight. Treatment and diet should only take place in consultation with your doctor.

Recommendations on exercise programmes or weight loss diets, on the basis of the values determined, should be provided by your doctor or certified health professional. Do not try to interpret the data yourself. SOEHNLE will not be liable for damage or loss caused by the Body Balance Adria or for third party claims.

- ▶ This product is intended for personal home use only. It is not designed for professional use in hospitals or medical facilities.
- ▶ Not suitable for people with electronic implants, such as pacemakers, etc.

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**Preparation for use**

- ① Insert batteries.  
The batteries should be removed if the scale is not going to be used for a while (see fig. ① below). Stored data are lost when the batteries are removed.
- ② Always place the scale on a level, firm surface for measurements.
- ③ Caution: slippery when wet.
- ④ Cleaning and care: clean only with a damp cloth. Never use solvents or abrasives. Do not immerse the scale in water.
- ⑤ Country-specific settings: the scale can be switched from kg/cm to lb/in or st/in (applies to the English-language versions only).

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






## Weighing operation

- ① Tap your foot lightly on the platform to start the scale.
- ② Wait a few seconds for **0.0** to appear in the display.
- ③ Step on to the scale and stand still until the weight symbol (♁) appears beside the display.
- ④ Body weight indication.
- ⑤ The scale switches off automatically when the person steps off.

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




## Individual programming

- ① Place scale on a table.
- ② Switch on scale
- ③ and start programming (↔). Select memory location (P1 ... P4) (1 key = minus, 2 key = plus) and confirm (↔)
- ④ Input your height (1 key = minus, 2 key = plus) and confirm (↔)
- ⑤ Input your age (1 key = minus, 2 key = plus) and confirm (↔)
- ⑥ Input your gender (1 key = ♀ male, 2 key = ♂ female), then press for 3 seconds to confirm (↔)
- ⑦ The scale is now programmed and will switch off automatically.

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









## Body fat/body water monitoring

**Always measure in bare feet.** Personal data must be programmed first (page 10).

- 1 Switch on scale and wait for **0.0** to appear in the display.
- 2 Select individual memory location (**key 1 ... 4**). **P1 ... P4** is shown in the display.
- 3 **0.0** and the gender symbol (**♂ / ♀**) appear. Scale is ready for use. Step on to the scale. Your weight is indicated, with a bar display while analysis is taking place (20-30 seconds). Stand still until analysis is complete.
- 4 After this, the body water (**≈**) and body fat (**↔**) content are displayed alternately.
- 5 Step off the scale. Display shows body weight.

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

## Interpretation of results

The Body Balance Adria determines your personal ideal weight and body water and body fat composition from your individual data and recommendations made by health experts, DGE (Deutsche Gesellschaft für Ernährung e.V., the German nutritional association) and the World Health Organisation (WHO).

The results of body water measurement are influenced by the proportion of body fat and muscle. If the proportion of body fat is high or the proportion of muscle is too low, the result of body water measurement will always be too low.

For consistent results we recommend that you always measure at the same time each day.

	Specific training recommended - water balance ok	Specific training recommended – make sure you drink plenty of fluid
	Everything OK: body water, body fat, muscle mass and weight are balanced	Body fat and weight are OK, but you are not drinking enough
	Change of diet - water balance ok	A change of diet is recommended – make sure you drink plenty of fluid

## The importance of water for your health

The human body is about 55-60% water – the amount varies according to age and gender.

Water performs a number of important roles in the body:

Water is a **constituent** of our cells. This means that all the cells in the body, whether in the skin, glands, muscles, brain or anywhere else, can only do their job if they have enough water.

Water acts as a **solvent** of substances important to the body.

Water is a **carrier** of nutrients, endogenous substances and products of metabolism.

Water is a **reactant** in chemical reactions involving or requiring water that take place continuously in each of the body's cells.

Water is a **coolant** and plays a part in regulating the temperature balance, particularly through perspiration.

The percentage of body water decreases with age. It can be:

- ▶ 75 to 80 per cent in newborn babies
- ▶ about 60 per cent in adult men of normal weight,
- ▶ 50 to 55 per cent in adult women of normal weight,
- ▶ about 50 per cent in older people.

Each of us loses a certain quantity of fluid every day.

The quantity can vary considerably and it needs to be regularly replaced.

Acute shortage of water can produce the following symptoms:

- ▶ Thirst, in the case of over 0.5 per cent of the body weight (0.3 to 0.4 litres of water)
- ▶ A decline in physical and mental performance is apparent on losing as little as 2 per cent of body water, for example through sweating or diarrhoea
- ▶ A 3 per cent drop causes a decline in production of saliva and urine, the mouth becomes dry and constipation may occur

- ▶ A greater loss of body water can lead to a rapid pulse, increased body temperature and confusion.

A healthy adult should normally drink about 2 litres of fluid. More fluid is needed in the case of vigorous expenditure of energy, heat, cold dry air, a high-salt diet, a high protein intake and when unwell, for example high temperature, vomiting, diarrhoea, etc. Always consult your doctor in the case of fluid loss due to illness.

### Health experts recommend the following percentages of body fat/body water:

Age	♂ Body fat				Body water	♀ Body fat				Body water		
	Excellent ++	Good +	Average •	High –		Exzellent ++	Good +	Average •	High –			
20–24	18.2	22.0	25.0	29.6	50–55	10.8	14.9	19.0	23.3	60		
25–29	18.9	22.1	25.4	29.8		12.8	16.5	20.3	24.3			
30–34	19.7	22.7	26.4	30.5		14.5	18.0	21.5	25.2			
35–39	21.0	24.0	27.7	31.5		16.1	19.3	22.6	26.1			
40–44	22.6	25.6	29.3	32.8		17.5	20.5	23.6	26.9			
45–49	24.3	27.3	30.9	34.1		18.6	21.5	24.5	27.6			
50–59	26.6	29.7	33.1	36.2		19.8	22.7	25.6	28.7			
> 60	27.4	30.7	34.0	37.3		50	20.2	23.7	26.2		29.3	50

Fluctuations can nevertheless occur (see page 19).

## Factors influencing the results

Analysis is based on the measurement of the body's electrical resistance. The water content in the lean mass and the proportion of body fat and muscle affect measurement results. Eating and drinking habits during the course of the day and individual lifestyle affect the water balance. This is noticeable due to fluctuations in body water and body fat percentage readings.

### Other factors can affect water balance too:

- ▶ After a bath, the body fat reading may be too low and the body water reading too high.
- ▶ After a meal, readings can be higher.
- ▶ Women may experience fluctuations due to the menstrual cycle.
- ▶ Loss of water caused by illness or after physical activity (sport). After exercise you should wait 6 to 8 hours before measuring your body fat percentage again.

It is important to always to measure yourself under the same conditions. **Varying or implausible results can occur in the case of:**

- Persons with a high temperature, symptoms of oedema or osteoporosis
- Persons on dialysis treatment
- Persons taking cardiovascular medicine
- Pregnant women
- Athletes who undertake more than 10 hours of intensive training per week and have a resting pulse rate of under 60/min
- Competitive athletes and body builders.

## Note

The scale has an energy-saving mode that switches it off automatically if no keys are pressed for approximately 1 minute. To continue programming, switch on the scale and keep pressing the confirmation key (↵) until data previously input appear. Then continue as described above.

### Technical data:

- ▶ Weighing capacity x graduation = max. 150 kg x 100 g
- ▶ Graduation: body fat: 0.1 %
- ▶ Graduation: body water: 0.1 %
- ▶ Height range: 100–250 cm
- ▶ Age range: 10–99 years
- ▶ 4-person memory for individual specific data storage
- ▶ Large LCD display
- ▶ Measuring current: 0.035 mA RMS, 50 kHz
- ▶ Battery requirements: 4 1.5 V Mignon AA (LR6) (batteries are included)

## Meaning of symbols

- ① Batteries spent.
- ② Overload: 150 kg or more
- ③ Implausible analysis result – check individual programming
- ④ Bad foot contact – clean scale or wash feet

## Warranty

Soehnle guarantees that all defects due to materials or manufacturing faults will be remedied by replacement or repair, free of charge, for a period of 3 years from the date of purchase. Please keep your purchase receipt and the guarantee card in a safe place. If you have any complaints, please return the scale to your dealer with the guarantee card and receipt.

### This device is screened in accordance with the applicable EC Directive 89/336/EEC.

Note: the display value may be affected by extreme electromagnetic influences, e.g. when a radio is operated in the immediate vicinity of the device. The product can be used for its intended purpose again when the interference disappears, but may have to be switched on again first.

### Correct battery disposal:

Batteries must not be disposed of as household waste. The law requires that you, as consumer, return the waste batteries either to public collection points in your town or to any outlet selling batteries of the same kind.

Note:

Batteries containing pollutant substances are marked as follows:  
 Pb = Battery contains lead  
 Cd = Battery contains cadmium  
 Hg = Battery contains mercury

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SOEHNLE-Waagen GmbH & Co. KG